

The **Debrisoft**[®] Difference




For removing debris and superficial slough from wounds or skin, there really is no better choice than Debrisoft.[®]

Patented Monofilament Fibre Technology™

- 18 million monofilament fibres
- Cut to a specific angle and length
- Soft and conformable – gentle on patients
- Angled tips reach uneven areas of the skin or wound bed

Unique mode of action

- **Lifts** up debris, superficial slough and exudate quickly
- **Binds** it within the pad
- **Removes** barriers to healing
- Leaves the area clear and promotes development of healthy tissue
- Will not damage any new granulation tissue and epithelial cells



RECOMMENDED BY  **NICE**.
DEBRISOFTING™ EFFECTIVE DEBRIDEMENT.
HEALTHY TISSUE.
THE DEBRISOFT[®] DIFFERENCE.

Debrisoft[®] Ordering information

Size (cm)	Ref	PIP code	NPC code	Bunzl	Credenhill	PECOS	Item number (Wales)	Units per pack
10x10	31222	358-1287	ELZ354	AL358-1287	31222	164271	ELC0281	5

How to use Debrisoft®

Step 1. Open the Debrisoft® single-use, sterile pack.

Step 2. Fully moisten the soft, fleecy side of Debrisoft® with 20–40ml (1–2 egg cups) of tap water or saline.*



Step 3. With gentle pressure, as if exfoliating, use the soft, fleecy side of the moistened Debrisoft® in a circular motion (on wounds) or in long, sweeping strokes (on the skin).



Step 4. Typically, a new piece of Debrisoft® is used for each separate wound/area of skin.* Dispose of the used Debrisoft® in normal clinical or household waste.

*Always refer to local guidelines.

Indications Debrisoft® is a rapid, highly effective, safe and easy method of debridement for acute and chronic wounds, peri-wound skin and skin conditions such as hyperkeratosis, in both adults and children. Debrisoft® can also be used to clear the wound bed in advance of the initial wound assessment or to aid accurate pressure ulcer categorisation.

Information and advice

- ✓ Prior to using Debrisoft®, a full holistic assessment should be completed.
- ✓ Wash off any emollients before using Debrisoft®.
- ✓ Always fully moisten Debrisoft® with 20–40ml of tap water or saline before use. (Always refer to local guidelines.)
- ✓ Always use the soft, fleecy side and not the knitted, reverse side.
- ✓ Debrisoft® is latex-free.
- ✓ Debrisoft® is for single use only.
- ✓ Where very stubborn slough or hard necrosis is present, treatment by autolytic debridement, using dressings (e.g. ActiformCool®) prior to treatment with Debrisoft® would be advised.
- ✓ Following the use of Debrisoft®, reassess the wound/skin before selecting an appropriate dressing/therapy.
- ✓ In some cases, it may be necessary to use Debrisoft® on more than one occasion.
- ✓ Larger areas may take longer to treat, depending on the amount of slough, debris or hyperkeratosis present.
- ✗ Debrisoft® must not be used as a wound dressing.
- ✗ Do not cut Debrisoft®.
- ✗ Debrisoft® should not be used if there is a known sensitivity to any components of the product (100% polyester).

ACTIVA®
HEALTHCARE
an **L&R** Company

Call our customer care line: **08450 606707**
International enquiries: **+44 1283 576800**
Or visit our website at: **www.debrisoft.co.uk**

Activa® and ActiFormCool® are registered trademarks of Activa Healthcare Ltd. Debrisoft® is a registered trademark of L&R. Monofilament Fibre Technology™ and Debrisofting™ are trademarks of L&R.

M1554 V1.4 September 2015

