

## ReadyWrap<sup>™</sup>

### Leg sizing chart and order form



New on Drug Tariff XXL calf,  
knee and thigh sizes  
(from August)

## ReadyWrap™

### Circumference points:

- cG** 2cm below gluteal fold
- cE<sup>1</sup>** 5cm above the knee
- cD** 2-3cm below popliteal fossa
- cC** Widest part of calf
- cB** 2cm above malleolus
- cY** Around instep and heel with foot dorsi-flexed
- cA** Around base of toes

### Lengths:

- IA-Y** Heel to base of toes
- ID-B** The length between point B (2cm above the malleolus) and point D (2cm below popliteal fossa)
- IG-E<sup>1</sup>** The length between point E<sub>1</sub> (5cm above the knee) and point G (2cm below gluteal fold)

**NB** **c** = circumference, **l** = length

all length measurements should be longitudinal and should follow body contours.

### How to measure the leg:

Ask the patient to stand, if possible.  
Mark and measure the leg as follows

### For the Thigh Garment:

**cG, cE<sup>1</sup>, IG-E<sup>1</sup>**

### For the Knee Garment:

**cE<sup>1</sup>, cD**

### For the Calf:

**cC, cB, ID-B**

### For the Foot Garment:

**cB, cY, IA-Y**

For further information on ReadyWrap garments and training please contact Activa Customer Care on **08450 606707** and visit our website [www.activahealthcare.co.uk/readywrap](http://www.activahealthcare.co.uk/readywrap) or contact your local Activa representative.

Delivery will be within 7 working days of receipt of a correctly completed order form; subject to circumstances beyond our control.

Please send the completed order form to Credenhill via:

Email: [sales@credenhill.co.uk](mailto:sales@credenhill.co.uk)

Fax: **0115 944 0437**

Post: **Credenhill Limited, 10 Cossall Industrial Estate,  
Ilkeston, Derbyshire, DE7 5UG.**

Phone: **0115 932 0144**



an  Company

Name of Pharmacy or Clinic: \_\_\_\_\_

Contact Name: \_\_\_\_\_

Account Number (if known): \_\_\_\_\_ Branch Number (if applicable): \_\_\_\_\_

Address: \_\_\_\_\_

Postcode: \_\_\_\_\_

Tel: \_\_\_\_\_ Email: \_\_\_\_\_


Practitioner Name: \_\_\_\_\_ Tel: \_\_\_\_\_

Order Ref. or Prescription No.: \_\_\_\_\_

Patient Name/Ref: \_\_\_\_\_

Date of Measurement: \_\_\_\_/\_\_\_\_/\_\_\_\_ Order Date: \_\_\_\_/\_\_\_\_/\_\_\_\_

**NB** All length measurements should follow body contour.  
All measurements to be in centimetres.



Write your measurement results on **bold line** adjacent to measuring point

Please tick where applicable to indicate your garment choice

Garment (Right Leg)	Length/Side	Colour	Size	Quantity	Code <small>transfer code from back of sheet</small>
Thigh	25cm 30cm <input type="checkbox"/> <input type="checkbox"/>	Beige Black <input type="checkbox"/> <input type="checkbox"/>			
Knee		Beige Black <input type="checkbox"/> <input type="checkbox"/>			
Calf	30cm 35cm <input type="checkbox"/> <input type="checkbox"/>	Beige Black <input type="checkbox"/> <input type="checkbox"/>			
Foot	Right <input type="checkbox"/>	Beige Black <input type="checkbox"/> <input type="checkbox"/>			

Garment (Left Leg)	Length/Side	Colour	Size	Quantity	Code <small>transfer code from back of sheet</small>
Thigh	25cm 30cm <input type="checkbox"/> <input type="checkbox"/>	Beige Black <input type="checkbox"/> <input type="checkbox"/>			
Knee		Beige Black <input type="checkbox"/> <input type="checkbox"/>			
Calf	30cm 35cm <input type="checkbox"/> <input type="checkbox"/>	Beige Black <input type="checkbox"/> <input type="checkbox"/>			
Foot	Left <input type="checkbox"/>	Beige Black <input type="checkbox"/> <input type="checkbox"/>			

Please tick where applicable to indicate your garment choice

<b>ReadyWrap™ Thigh (IG-E<sup>1</sup>)</b> (2 x lengths)	<b>Small</b> Groin Circumference (cG) 48-63cm Above Knee Circumference (cE <sup>1</sup> ) 38-48cm	<b>Medium</b> Groin Circumference (cG) 58-73cm Above Knee Circumference (cE <sup>1</sup> ) 45-58cm	<b>Large</b> Groin Circumference (cG) 68-83cm Above Knee Circumference (cE <sup>1</sup> ) 52-68cm	<b>Extra Large</b> Groin Circumference (cG) 78-93cm Above Knee Circumference (cE <sup>1</sup> ) 59-78cm	<b>Extra Extra Large</b> Groin Circumference (cG) 88-103cm Above Knee Circumference (cE <sup>1</sup> ) 66-88cm
	<input type="checkbox"/> 104478 <input type="checkbox"/> 104479 <input type="checkbox"/> 104480 <input type="checkbox"/> 104481	<input type="checkbox"/> 104474 <input type="checkbox"/> 104475 <input type="checkbox"/> 104476 <input type="checkbox"/> 104477	<input type="checkbox"/> 104470 <input type="checkbox"/> 104471 <input type="checkbox"/> 104472 <input type="checkbox"/> 104473	<input type="checkbox"/> 104482 <input type="checkbox"/> 104483 <input type="checkbox"/> 104484 <input type="checkbox"/> 104485	<input type="checkbox"/> 104790 <input type="checkbox"/> 104791 <input type="checkbox"/> 104792 <input type="checkbox"/> 104793
<b>ReadyWrap™ Knee</b> (one length - 20cm)	<b>Small</b> Above Knee Circumference (cE <sup>1</sup> ) 38-48cm Below Knee Circumference (cD) 28-38cm	<b>Medium</b> Above Knee Circumference (cE <sup>1</sup> ) 45-58cm Below Knee Circumference (cD) 33-43cm	<b>Large</b> Above Knee Circumference (cE <sup>1</sup> ) 52-68cm Below Knee Circumference (cD) 38-48cm	<b>Extra Large</b> Above Knee Circumference (cE <sup>1</sup> ) 59-78cm Below Knee Circumference (cD) 43-58cm	<b>Extra Extra Large</b> Above Knee Circumference (cE <sup>1</sup> ) 66-88cm Below Knee Circumference (cD) 53-68cm
	<input type="checkbox"/> 104464 <input type="checkbox"/> 104465	<input type="checkbox"/> 104462 <input type="checkbox"/> 104463	<input type="checkbox"/> 104460 <input type="checkbox"/> 104461	<input type="checkbox"/> 104466 <input type="checkbox"/> 104467	<input type="checkbox"/> 104788 <input type="checkbox"/> 104789
<b>ReadyWrap™ Calf (ID-B)</b> (2 x lengths)	<b>Small</b> Calf Circumference (cC) 28-38cm Ankle Circumference (cB) 18-23cm	<b>Medium</b> Calf Circumference (cC) 33-43cm Ankle Circumference (cB) 20-28cm	<b>Large</b> Calf Circumference (cC) 38-48cm Ankle Circumference (cB) 25-33cm	<b>Extra Large</b> Calf Circumference (cC) 43-58cm Ankle Circumference (cB) 30-38cm	<b>Extra Extra Large</b> Calf Circumference (cC) 53-68cm Ankle Circumference (cB) 35-43cm
	<input type="checkbox"/> 104450 <input type="checkbox"/> 104451 <input type="checkbox"/> 104452 <input type="checkbox"/> 104453	<input type="checkbox"/> 104446 <input type="checkbox"/> 104447 <input type="checkbox"/> 104448 <input type="checkbox"/> 104449	<input type="checkbox"/> 104442 <input type="checkbox"/> 104443 <input type="checkbox"/> 104444 <input type="checkbox"/> 104445	<input type="checkbox"/> 104454 <input type="checkbox"/> 104455 <input type="checkbox"/> 104456 <input type="checkbox"/> 104457	<input type="checkbox"/> 104784 <input type="checkbox"/> 104785 <input type="checkbox"/> 104786 <input type="checkbox"/> 104787
<b>ReadyWrap™ Foot (Side)</b>	<b>Small</b> Foot Length (IA-Y) up to 21cm Ankle Circumference (cB) 18-23cm Circumference (cY) 28-33cm	<b>Medium</b> Foot Length (IA-Y) up to 23cm Ankle Circumference (cB) 20-28cm Circumference (cY) 30-38cm	<b>Large</b> Foot Length (IA-Y) up to 26cm Ankle Circumference (cB) 25-33cm Circumference (cY) 35-43cm	<b>Extra Large</b> Foot Length (IA-Y) up to 29cm Ankle Circumference (cB) 30-38cm Circumference (cY) 40-48cm	
	<b>Right</b> <input type="checkbox"/> Beige <b>Right</b> <input type="checkbox"/> Black <b>Left</b> <input type="checkbox"/> Beige <b>Left</b> <input type="checkbox"/> Black	<input type="checkbox"/> 104436 <input type="checkbox"/> 104437 <input type="checkbox"/> 104428 <input type="checkbox"/> 104429	<input type="checkbox"/> 104434 <input type="checkbox"/> 104435 <input type="checkbox"/> 104426 <input type="checkbox"/> 104427	<input type="checkbox"/> 104432 <input type="checkbox"/> 104433 <input type="checkbox"/> 104424 <input type="checkbox"/> 104425	<input type="checkbox"/> 104438 <input type="checkbox"/> 104439 <input type="checkbox"/> 104430 <input type="checkbox"/> 104431

Delivery will be within 7 working days of receipt of a correctly completed order form; subject to circumstances beyond our control.

Please send the completed order form to Credenhill via:

Email: [sales@credenhill.co.uk](mailto:sales@credenhill.co.uk)

Fax: 0115 944 0437

Post: Credenhill Limited, 10 Cossall Industrial Estate, Ilkeston, Derbyshire, DE7 5UG.

Phone: 0115 932 0144