

## measuring & fitting information

**A:** 5cm below the groin

**B:** Widest part of the calf

**C:** Just above the malleolus (ankle bone)

**D:** From back of the heel to the longest toe

| Activa® British Standard hosiery measuring guide                   |             |             |             |             |             |
|--|-------------|-------------|-------------|-------------|-------------|
| a size indicator can be found by a coloured thread in the stocking |             |             |             |             |             |
| Size   | S           | M           | L           | XL          | XXL         |
| A  | 40.5-53.0cm | 43.0-55.0cm | 48.0-60.5cm | 53.0-66.0cm | N/A         |
| B  | 30.5-37.5cm | 33.0-40.0cm | 35.5-43.0cm | 38.0-46.0cm | 40.0-49.0cm |
| C  | 19.5-25.5cm | 21.5-27.5cm | 22.5-29.5cm | 23.0-32.0cm | 23.5-34.0cm |
| D  | 20.5-24.0cm | 23.0-26.0cm | 25.5-29.5cm | 26.5-32.5cm | 26.5-32.5cm |

