

Activa® British Standard Hosiery - Measure and Sizing Form

Patient name _____


Date _____

Clinician _____

Clinic or Pharmacy Address:

(Please specify the address where the patient will collect their hosiery)

Step 1 - Measure

	Left leg	Right leg
 A: 5cm below groin		
B: Widest part of the calf		
C: Just above the malleolus (ankle bone)		
D: From back of heel to the longest toe		

COMPRESSION REMINDER SERVICE

Encourage the wearer to **scan this QR code** and register for the **L&R Compression Reminder Service**.



They will receive timely reminders about scheduling their hosiery renewal appointment.

	Size requested
	Number of packs

Is the ActiGlide® Hosiery Applicator required?

(PIP Code: 284-5196) **Yes** **No**

Step 2 - Select size

	Red	Blue	Green	Yellow	Purple
Size	Small	Medium	Large	Extra large	Extra Extra large*
A	40.5 – 53.0cm	43.0 – 55.0cm	48.0 – 60.5cm	53.0 – 66.0cm	N/A
B	30.5 – 37.5cm	33.0 – 40.0 cm	35.5 – 43.0cm	38.0 – 46.0cm	40.0 – 49.0cm
C	19.5 – 25.5cm	21.5 – 27.5cm	22.5 – 29.5cm	23.0 – 32.0cm	23.5 – 34.0cm
D	20.5 – 24.0cm	23.0 – 26.0cm	25.5 – 29.5cm	26.5 – 32.5cm	26.5 – 32.5cm

* For availability of XXL garments, please see PIP codes on reverse.

Step 3 - Select garment

Please tick a box to indicate garment required	Colour	10mmHg	14-17mmHg Class 1	18-24mmHg Class 2	25-35mmHg Class 3	40mmHg
Compression Liner (Closed Toe)†	White					
	Sand					
	Black					
Compression Liner (Open Toe)	Sand					
Below Knee Stockings (Closed Toe)	Sand					
	Honey					
	Black					
Below Knee Stockings (Open Toe)	Sand					
Thigh Length Stockings (Closed Toe)	Sand					
	Honey					
	Black					
Thigh Length Stockings (Open Toe)	Sand					
Unisex Sock (Closed Toe)	Black					
	Brown					
Unisex Patterned Sock (Closed Toe)	Black					
Leg Ulcer Hosiery Kit (1 x below knee open toe stocking plus 2 compression liners)	Sand Stocking/White Liners					
	Black Stocking/Sand Liners					

† Search under 'Activa Leg Ulcer Hosiery Kit' to find Activa Compression Liners on your ordering system.

Activa® British Standard Hosiery Range

Product & Style Colour of size indicator thread in stocking	Colour	Small Red	Medium Blue	Large Green	Extra Large Yellow	Extra Extra Large Purple
Liner Pack - 10mmHg Provides light compression for maintenance therapy, palliative care. They may be used with other compression stockings to increase compression on the limb, or to build the compression level in layers if the patient struggles to apply a higher class of hosiery.						
Stocking Liner (3 Singles)	White	311-8767	311-8775	311-8783	311-8791	311-8809
Stocking Liner (3 Singles)	Sand	326-9750	326-9768	326-9776	326-9784	326-9792
Stocking Liner (3 Singles)	Black	387-0805	387-0813	387-0821	387-0839	387-0847
Stocking Liner Open Toe (3 Singles)	Sand	340-3557	340-3565	340-3573	340-3581	340-3599
Class 1 - 14-17mmHg Provides effective compression in the treatment of superficial and early venous disease progression and varicose veins, including those during pregnancy.						
Thigh Length	Sand	259-0305	259-0339	259-0362	259-0396	-
Thigh Length	Honey	259-0313	259-0347	259-0370	259-0404	-
Thigh Length	Black	259-0321	259-0354	259-0388	259-0412	-
Below Knee	Sand	259-0420	259-0446	259-0461	259-0487	-
Below Knee	Honey	259-0438	259-0453	259-0479	259-0495	-
Below Knee	Black	278-2373	278-2381	278-2399	278-2407	-
Below Knee Open Toe	Sand	268-1500	268-1518	268-1526	268-1534	317-2723
Unisex Sock	Black	278-2415	278-2423	278-2431	278-2449	-
Unisex Patterned Sock	Black	356-9993	357-0009	357-0017	357-0025	-
Class 2 - 18-24mmHg Provides compression for medium severity varicose veins or prevention of venous leg ulcers.						
Thigh Length	Sand	259-0503	259-0537	259-0560	259-0594	-
Thigh Length	Honey	259-0511	259-0545	259-0578	259-0602	-
Thigh Length	Black	259-0529	259-0552	259-0586	259-0610	-
Thigh Length Open Toe	Sand	259-0859	259-0867	259-0875	259-0883	-
Below Knee	Sand	259-0628	259-0685	259-0743	259-0784	-
Below Knee	Honey	259-0651	259-0701	259-0750	259-0792	-
Below Knee	Black	259-0669	259-0735	259-0776	259-0800	-
Below Knee Open Toe	Sand	259-0818	259-0826	259-0834	259-0842	317-2707
Unisex Sock	Black	259-4463	259-0891	259-0917	259-0933	-
Unisex Sock	Brown	259-4471	259-0909	259-0925	259-0941	-
Unisex Patterned Sock	Black	357-0033	357-0041	357-0058	357-0066	-
Class 3 - 25-35mmHg Provides effective management of chronic venous insufficiency and eases away the pain and discomfort of severe varicose veins and chronic leg ulcers. The hosiery can also be used in the prevention of leg ulcer recurrence.						
Thigh Length Open Toe	Sand	259-1287	259-1295	259-1303	259-1311	-
Below Knee Open Toe	Sand	259-1246	259-1253	259-1261	259-1279	310-2977
40mmHg Activa® Leg Ulcer Hosiery Kit Designed to be a treatment and maintenance solution in the management of uncomplicated venous leg ulcers.						
Sand (1 Sand Class 3 Stocking + 2 White Liners)	-	311-9021	311-9047	311-9054	311-9070	311-9088
Black (1 Black Class 3 Stocking + 2 Sand Liners)	-	366-5593	366-5601	366-5619	366-5627	366-5635
ActiGlide® Hosiery Applicator						
ActiGlide® Hosiery Applicator	One Size Fits All	284-5196				



For further information please call **08450 606 707**
to speak directly with our Customer Services Team
or visit our website **www.Lohmann-Rauscher.co.uk**

Activa®, ActiLymph® and ActiGlide® are registered to L&R Medical UK Ltd.



M002 V4.8