

ReadyWrap® Measuring the full lower limb only

For the below knee measuring form if oedema management below upper limb


Scan the QR code to access the form ▶



Measuring form for full lower limb - Complete one form per leg

Each garment comes with 1 x ReadyWrap® Liner

How to measure the leg: Measure the leg whilst patient is standing, if possible.



Write your measurement results in the box adjacent to measuring point

Please note

c = circumference, L = length

All length measurements should be longitudinal and should follow body contours.

All measurements to be in centimetres.

Patients Name: _____ Clinicians Name: _____

Date of measuring: _____ Contact Number: _____

Delivery: to orderer directly to patient

| ReadyWrap® Thigh | | | | | |
|---|-------------------|----------|----------|-------------------|-----------|
| Circumference c / Length L | Small | Medium | Large | XL | XXL |
| cG 2cm below gluteal fold | 48-63 cm | 58-73 cm | 68-83 cm | 78-93 cm | 88-103 cm |
| cE¹ 5cm above the knee | 38-48cm | 45-58cm | 52-68cm | 59-78cm | 66-88cm |
| LG-E¹ the length between point cE¹ (5cm above the knee) and point cG (2cm below the gluteal fold) | Up to 25cm – 25cm | | | Up to 30cm – 30cm | |

| ReadyWrap® Knee | | | | | |
|--------------------------------|---------|---------|---------|---------|---------|
| Circumference c / Length L | Small | Medium | Large | XL | XXL |
| cD 2-3cm below popliteal fossa | 28-38cm | 33-43cm | 38-48cm | 43-58cm | 53-68cm |
| cE¹ 5cm above the knee | 38-48cm | 45-58cm | 52-68cm | 59-78cm | 66-88cm |

| ReadyWrap® Calf | | | | | |
|---|-------------------|---------|---------|-------------------|---------|
| Circumference c / Length L | Small | Medium | Large | XL | XXL |
| cB 2cm above malleolus | 18-23cm | 20-28cm | 25-33cm | 30-38cm | 35-43cm |
| cC widest part of calf | 28-38cm | 33-43cm | 38-48cm | 43-58cm | 53-68cm |
| LD-B the length between point B (2cm above the malleolus) and 2cm below popliteal fossa | Up to 30cm – 30cm | | | Up to 35cm – 35cm | |

| ReadyWrap® Foot CT | | | | | |
|--|----------------------|-------------|-----------|----------------|-----------|
| Circumference c / Length L | Small | Medium | Large | XL | XXL |
| cA around base of the toes | 18-20.4cm | 20.5-23.5cm | 23.6-26cm | 26.1-28.5cm | 28.6-31cm |
| cA¹ widest part around the dorsum of the foot | 20-25cm | 24-29cm | 28-33cm | 32-37cm | 36-41cm |
| cB 2cm above malleolus | 18-23cm | 20-28cm | 25-33cm | 30-38cm | 35-43cm |
| LA-Y from the base of the little toe to the back of the heel | Up to 19cm – Regular | | | 19-22cm – Long | |

| ReadyWrap® Toe | | | | | |
|----------------------------|-----------|-------------|-----------|-------------|-----------|
| Circumference c / Length L | Small | Medium | Large | XL | XXL |
| cA around base of the toes | 18-20.4cm | 20.5-23.5cm | 23.6-26cm | 26.1-28.5cm | 28.6-31cm |

| | | |
|--------------------------|-------|-------|
| Colour | Black | Beige |
| Left or Right Leg | Left | Right |

| | | |
|--|--------------------------|---------------------|
| ReadyWrap® Extender Strap | Black | Beige |
| Additional ReadyWrap® Liners | Small – Large cB 18-33cm | XL – XXL cB 30-43cm |
| Sold in pairs, available in black only | | |

Correct ordering description example:

ReadyWrap® Foot CT Small Regular Black Left x1 / ReadyWrap® Calf Small 30cm (x1) Black

Please tick where applicable to indicate your garment choice:

| | Small | Medium | Large | XL | XXL |
|--|--------|--------|--------|--------|--------|
| ReadyWrap® Toe | | | | | |
| Right ■ Beige | 104416 | 104414 | 104412 | 104418 | 104420 |
| Right ■ Black | 104417 | 104415 | 104413 | 104419 | 104421 |
| Left ■ Beige | 104406 | 104404 | 104402 | 104408 | 104410 |
| Left ■ Black | 104407 | 104405 | 104403 | 104409 | 104411 |
| ReadyWrap® Foot CT | | | | | |
| Right ■ Beige Reg up to 19cm | 141727 | 141731 | 141735 | 141739 | 141743 |
| Right ■ Black Reg up to 19cm | 141729 | 141733 | 141737 | 141741 | 141745 |
| Left ■ Beige Reg up to 19cm | 141726 | 141730 | 141734 | 141738 | 141742 |
| Left ■ Black Reg up to 19cm | 141728 | 141732 | 141736 | 141740 | 141744 |
| Right ■ Beige Long 19-22cm | 141747 | 141751 | 141755 | 141759 | 141763 |
| Right ■ Black Long 19-22cm | 141749 | 141753 | 141757 | 141761 | 141765 |
| Left ■ Beige Long 19-22cm | 141746 | 141750 | 141754 | 141758 | 141762 |
| Left ■ Black Long 19-22cm | 141748 | 141752 | 141756 | 141760 | 141764 |
| ReadyWrap® Calf | | | | | |
| ■ Beige up to 30 cm | 104450 | 104446 | 104442 | 104454 | 104784 |
| ■ Black up to 30 cm | 104451 | 104447 | 104443 | 104455 | 104785 |
| ■ Beige up to 35 cm | 104452 | 104448 | 104444 | 104456 | 104786 |
| ■ Black up to 35 cm | 104453 | 104449 | 104445 | 104457 | 104787 |
| ReadyWrap® Knee | | | | | |
| ■ Beige | 104464 | 104462 | 104460 | 104466 | 104788 |
| ■ Black | 104465 | 104463 | 104461 | 104467 | 104789 |
| ReadyWrap® Thigh | | | | | |
| ■ Beige up to 25 cm | 104478 | 104474 | 104470 | 104482 | 104790 |
| ■ Black up to 25 cm | 104479 | 104475 | 104471 | 104483 | 104791 |
| ■ Beige up to 30 cm | 104480 | 104476 | 104472 | 104484 | 104792 |
| ■ Black up to 30 cm | 104481 | 104477 | 104473 | 104485 | 104793 |

ReadyWrap® Liner Sold in pairs, available in black only

| | Small - Large | X-Large - XX-Large |
|--|---------------|--------------------|
| Below Knee ■ Black | 149032 | 149033 |
| Thigh High ■ Black | 149034 | 149035 |

ReadyWrap® Extender Strap One size fits all

| | |
|--|--------|
| ■ Beige | 104740 |
| ■ Black | 104741 |

Scan to download the **Compression Selector App**



Scan here to access **EPS and DAC** ordering support



Pharmacy Support
Scan here to access free pharmacy resources and information



Scan here to access our **FREE learning modules and videos**



Scan here to contact our **Customer Services Team**

