

# ActiLymph®

## Armsleeve Patient Record

This is **NOT** a made-to-measure form

To Doctor: \_\_\_\_\_

To the Pharmacist: \_\_\_\_\_

### Request for the ActiLymph® armsleeve

Patient Name: \_\_\_\_\_ Date: \_\_\_\_\_

Nurse/Clinic: \_\_\_\_\_ Name/Stamp: \_\_\_\_\_

#### European Class:

European Class 1 (18-21mmHg)

European Class 2 (23-32mmHg)

#### Style of armsleeve:

Armsleeve (CG)

Combined with handpiece (AG)

TopBand

#### Quantity:

1


2

3

4

N.B. Normal wear of armsleeve is 3 months (or 100 washes)

ActiGlide® hosiery applicator (PIP Code 284-5196)



**G:** The armpit crease

**E:** The elbow

**C:** The wrist

**A:** The widest part of the hand

#### ActiLymph® measuring guide for armsleeves

SIZE	S	M	L	XL
<b>G</b>	26.5-30.5cm	30.5-35cm	35-40cm	40-45cm
<b>E</b>	23.5-27cm	27-31cm	31-35.5cm	35.5-41cm
<b>C</b>	15-17cm	17-19.5cm	19.5-22.5cm	22.5-26cm
<b>A</b>	17.5-19.5cm	19.5-22cm	22-25cm	23-26cm

**Standard length**

**C-G**

<44cm

**Long length\***

**C-G**

>44cm

Before applying armsleeves, please read patient information leaflet enclosed within the box.

Left arm

Right arm

\_\_\_\_\_ **G** \_\_\_\_\_

\_\_\_\_\_ **E** \_\_\_\_\_

\_\_\_\_\_ **C** \_\_\_\_\_

\_\_\_\_\_ **A** \_\_\_\_\_

\_\_\_\_\_ **Size Requested** \_\_\_\_\_

\_\_\_\_\_ **Length Requested** \_\_\_\_\_

\* Only available with TopBand

**ActiLymph® European Class armsleeve range**
**Sizes & PIP Codes**

Product Description (Colour - Sand)	Small	Medium	Large	ExLarge
Class 1 Standard (CG)	317-3150	317-3309	317-3317	337-6928
Class 1 Standard (CG) TopBand	317-3325	317-3333	317-3341	336-8917
Class 1 Standard (CG) TopBand Long	342-8794	342-8802	342-8810	342-8828
Class 1 Combined Standard (AG)	317-3358	317-3366	317-3374	-
Class 1 Combined Standard (AG) TopBand	317-3382	317-3390	317-3408	-
Class 1 Combined Standard (AG) TopBand Long	342-8836	342-8844	342-8851	-
Class 2 Standard (CG)	317-3416	317-3424	317-3432	336-9345
Class 2 Standard (CG) TopBand	317-3440	317-3457	317-3465	336-9352
Class 2 Standard (CG) TopBand Long	342-8885	342-8893	342-8901	342-8919
Class 2 Combined Standard (AG)	317-3473	317-3481	317-3499	-
Class 2 Combined Standard (AG) TopBand	317-3507	317-3515	317-3523	-
Class 2 Combined Standard (AG) TopBand Long	342-8927	342-8935	342-8943	-



**ActiGlide® Hosiery Applicator**  
(PIP Code: 284-5196)

Call our customer care line: **08450 606 707**  
 International enquiries: **+44 1283 576 800**  
 or visit our website at: **[www.Lohmann-Rauscher.co.uk](http://www.Lohmann-Rauscher.co.uk)**