

ReadyWrap®

Measuring the upper limb.

How to measure the arm:

Measure the arm in a supported and relaxed position with the elbow slightly flexed.

The measurement points refer to the following anatomical landmarks:

Circumference points:

- cA** widest width of the Palm, just below the metacarpal phalangeal joints
- cC** narrowest point of the Wrist
- cE** around the Elbow joint and into the crease of the elbow
- cG** top of the Arm, just below the axillary crease


Lengths:

- IA-C** the length between point A (widest width of the palm, just below the metacarpal phalangeal joints) and point C (narrowest point of the wrist)
- IC-G** the length between point C (narrowest point of the wrist) and point G (top of the arm, just below the axillary crease)

NB c = circumference, l = length

All length measurements should follow body contour.

All measurements to be in centimetres.



Write your measurement results in the **box** adjacent to measuring point

Delivery: to orderer directly to patient

Name and address of recipient:

Patient Reference _____

Name/Last Name _____

Clinicians Name _____

Street _____

Clinicians Work Base _____

Postal Code/Place _____

Phone _____

Phone _____

Date of Measurement _____

Order Date _____

Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	Extra Large
ReadyWrap® Armsleeve				
(I-C-G) (3 x lengths)	Axilla Circumference (cG) 22-29 cm Elbow Circumference (cE) 19-26 cm Wrist Circumference (cC) 13.5-17 cm	Axilla Circumference (cG) 26-36 cm Elbow Circumference (cE) 24-33 cm Wrist Circumference (cC) 15.5-20 cm	Axilla Circumference (cG) 31-43 cm Elbow Circumference (cE) 29-40 cm Wrist Circumference (cC) 17.5-23 cm	Axilla Circumference (cG) 36-50 cm Elbow Circumference (cE) 30-43 cm Wrist Circumference (cC) 19.5-26 cm
Short				
■ Beige up to 43 cm	104514	104508	104502	104520
■ Black up to 43 cm	104515	104509	104503	104521
Average				
■ Beige up to 46 cm	104516	104510	104504	104522
■ Black up to 46 cm	104517	104511	104505	104523
Tall				
■ Beige up to 49.5 cm	104518	104512	104506	104524
■ Black up to 49.5 cm	104519	104513	104507	104525
ReadyWrap® Gauntlet				
	Length (IA-C) up to 12.5 cm Wrist Circumference (cC) 13.5-17 cm Palm Circumference (cA) 18-20 cm	Length (IA-C) up to 13.5 cm Wrist Circumference (cC) 15.5-20 cm Palm Circumference (cA) 20-22 cm	Length (IA-C) up to 14.5 cm Wrist Circumference (cC) 17.5-23 cm Palm Circumference (cA) 22-24 cm	Length (IA-C) up to 15.5 cm Wrist Circumference (cC) 19.5-26 cm Palm Circumference (cA) 24-26 cm
Right ■ Beige	104498	104496	104494	104500
Right ■ Black	104499	104497	104495	104501
Left ■ Beige	104490	104488	104486	104492
Left ■ Black	104491	104489	104487	104493

Please mark (tick or cross) where applicable to indicate the garment choice and then enter the sizing, quantity and reference code.

Garment (Armsleeve)	Length/Side	Colour	Size	Quantity	Code transfer code from back of sheet
Arm	43 cm	Beige Black			
	46 cm				
	49 cm				
Gauntlet	Right	Beige Black			
	Left				