


ReadyWrap® with Foot CT

Measuring the full lower limb.

How to measure the leg: Ask the patient to stand, if possible. Mark and measure the leg as follows



Right Leg

Left Leg

If patient bilateral, complete 2 forms

Write your measurement results in the box adjacent to measuring point

Toe

cA around base of the toes

Foot

cA around the base of the toes

cA' widest part around the dorsum of the foot

cB 2cm above malleolus

IA-Y from the base of the little toe to the back of the heel

Calf

cB 2cm above malleolus

cC widest part of calf

ID-B the length between point B (2cm above the malleolus) and point D (2cm below popliteal fossa)

Knee

cD 2-3cm below popliteal fossa

cE' 5cm above the knee

Thigh

cE' 5cm above the knee

cG 2cm below gluteal fold

IG-E' the length between point cE1 (5cm above the knee) and point cG (2cm below the gluteal fold)

Please note

c = circumference, l = length

All length measurements should be longitudinal and should follow body contours.

All measurements to be in centimetres.

Delivery: to orderer directly to patient Patient Reference _____

Clinicians Name _____ Clinicians Work Base _____

Phone _____ Date of Measurement _____

Name and address of recipient

Name / Last Name _____

Street _____

Postal Code _____

Phone _____ Order Date _____

Please mark (tick or cross) where applicable to indicate the garment choice and then enter the sizing, quantity and reference code.

Garment	Length	Foot	Colour	Size	Quantity	Code
Thigh	<input type="checkbox"/> 25 cm <input type="checkbox"/> 30 cm		<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Knee			<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Calf	<input type="checkbox"/> 30 cm <input type="checkbox"/> 35 cm		<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Foot CT	<input type="checkbox"/> Regular <input type="checkbox"/> Long	<input type="checkbox"/> Left <input type="checkbox"/> Right	<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Toe		<input type="checkbox"/> Left <input type="checkbox"/> Right	<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Extender Strap			<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Liner	<input type="checkbox"/> Below Knee <input type="checkbox"/> Thigh High		<input type="checkbox"/> Black			

Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	Extra Large	Extra Extra Large
ReadyWrap® Foot CT NEW					
(I-A-Y) (2 x lengths) (Side)	Ankle Circumference (cB) 18-23 cm Circumference (cA) 18-20.4 cm Circumference (cA1) 20-25 cm	Ankle Circumference (cB) 20-28 cm Circumference (cA) 20.5-23.5 cm Circumference (cA1) 24-29 cm	Ankle Circumference (cB) 25-33 cm Circumference (cA) 23.6-26 cm Circumference (cA1) 28-33 cm	Ankle Circumference (cB) 30-38 cm Circumference (cA) 26.1-28.5 cm Circumference (cA1) 32-37 cm	Ankle Circumference (cB) 35-43 cm Circumference (cA) 28.6-31 cm Circumference (cA1) 36-41 cm
Right <input type="checkbox"/> Beige Reg up to 19cm Right <input type="checkbox"/> Black Reg up to 19cm Left <input type="checkbox"/> Beige Reg up to 19cm Left <input type="checkbox"/> Black Reg up to 19cm Right <input type="checkbox"/> Beige Long 19-22cm Right <input type="checkbox"/> Black Long 19-22cm Left <input type="checkbox"/> Beige Long 19-22cm Left <input type="checkbox"/> Black Long 19-22cm	<input type="checkbox"/> 141727 <input type="checkbox"/> 141729 <input type="checkbox"/> 141726 <input type="checkbox"/> 141728 <input type="checkbox"/> 141747 <input type="checkbox"/> 141749 <input type="checkbox"/> 141746 <input type="checkbox"/> 141748	<input type="checkbox"/> 141731 <input type="checkbox"/> 141733 <input type="checkbox"/> 141730 <input type="checkbox"/> 141732 <input type="checkbox"/> 141751 <input type="checkbox"/> 141753 <input type="checkbox"/> 141750 <input type="checkbox"/> 141752	<input type="checkbox"/> 141735 <input type="checkbox"/> 141737 <input type="checkbox"/> 141734 <input type="checkbox"/> 141736 <input type="checkbox"/> 141755 <input type="checkbox"/> 141757 <input type="checkbox"/> 141754 <input type="checkbox"/> 141756	<input type="checkbox"/> 141739 <input type="checkbox"/> 141741 <input type="checkbox"/> 141738 <input type="checkbox"/> 141740 <input type="checkbox"/> 141759 <input type="checkbox"/> 141761 <input type="checkbox"/> 141758 <input type="checkbox"/> 141760	<input type="checkbox"/> 141743 <input type="checkbox"/> 141745 <input type="checkbox"/> 141742 <input type="checkbox"/> 141744 <input type="checkbox"/> 141763 <input type="checkbox"/> 141765 <input type="checkbox"/> 141762 <input type="checkbox"/> 141764
ReadyWrap® Thigh					
(I-G-E1) (2 x lengths)	Groin Circumference (cG) 48-63 cm Above Knee Circumference (cE) 38-48 cm	Groin Circumference (cG) 58-73 cm Above Knee Circumference (cE) 45-58 cm	Groin Circumference (cG) 68-83 cm Above Knee Circumference (cE) 52-68 cm	Groin Circumference (cG) 78-93 cm Above Knee Circumference (cE) 59-78 cm	Groin Circumference (cG) 88-103 cm Above Knee Circumference (cE) 66-88 cm
<input type="checkbox"/> Beige up to 25cm <input type="checkbox"/> Black up to 25cm <input type="checkbox"/> Beige up to 30cm <input type="checkbox"/> Black up to 30cm	<input type="checkbox"/> 104478 <input type="checkbox"/> 104479 <input type="checkbox"/> 104480 <input type="checkbox"/> 104481	<input type="checkbox"/> 104474 <input type="checkbox"/> 104475 <input type="checkbox"/> 104476 <input type="checkbox"/> 104477	<input type="checkbox"/> 104470 <input type="checkbox"/> 104471 <input type="checkbox"/> 104472 <input type="checkbox"/> 104473	<input type="checkbox"/> 104482 <input type="checkbox"/> 104483 <input type="checkbox"/> 104484 <input type="checkbox"/> 104485	<input type="checkbox"/> 104790 <input type="checkbox"/> 104791 <input type="checkbox"/> 104792 <input type="checkbox"/> 104793
ReadyWrap® Knee					
(one length - 20 cm)	Above Knee Circumference (cE) 38-48 cm Below Knee Circumference (cD) 28-38 cm	Above Knee Circumference (cE) 45-58 cm Below Knee Circumference (cD) 33-43 cm	Above Knee Circumference (cE) 52-68 cm Below Knee Circumference (cD) 38-48 cm	Above Knee Circumference (cE) 59-78 cm Below Knee Circumference (cD) 43-58 cm	Above Knee Circumference (cE) 66-88 cm Below Knee Circumference (cD) 53-68 cm
<input type="checkbox"/> Beige <input type="checkbox"/> Black	<input type="checkbox"/> 104464 <input type="checkbox"/> 104465	<input type="checkbox"/> 104462 <input type="checkbox"/> 104463	<input type="checkbox"/> 104460 <input type="checkbox"/> 104461	<input type="checkbox"/> 104466 <input type="checkbox"/> 104467	<input type="checkbox"/> 104788 <input type="checkbox"/> 104789
ReadyWrap® Calf					
(I-D-B) (2 x lengths)	Calf Circumference (cC) 28-38 cm Ankle Circumference (cB) 18-23 cm	Calf Circumference (cC) 33-43 cm Ankle Circumference (cB) 20-28 cm	Calf Circumference (cC) 38-48 cm Ankle Circumference (cB) 25-33 cm	Calf Circumference (cC) 43-58 cm Ankle Circumference (cB) 30-38 cm	Calf Circumference (cC) 53-68 cm Ankle Circumference (cB) 35-43 cm
<input type="checkbox"/> Beige up to 30 cm <input type="checkbox"/> Black up to 30 cm <input type="checkbox"/> Beige up to 35 cm <input type="checkbox"/> Black up to 35 cm	<input type="checkbox"/> 104450 <input type="checkbox"/> 104451 <input type="checkbox"/> 104452 <input type="checkbox"/> 104453	<input type="checkbox"/> 104446 <input type="checkbox"/> 104447 <input type="checkbox"/> 104448 <input type="checkbox"/> 104449	<input type="checkbox"/> 104442 <input type="checkbox"/> 104443 <input type="checkbox"/> 104444 <input type="checkbox"/> 104445	<input type="checkbox"/> 104454 <input type="checkbox"/> 104455 <input type="checkbox"/> 104456 <input type="checkbox"/> 104457	<input type="checkbox"/> 104784 <input type="checkbox"/> 104785 <input type="checkbox"/> 104786 <input type="checkbox"/> 104787
ReadyWrap® Toes					
(cA) (Sides)	(cA) 18.0-20.4 cm Circumference	(cA) 20.5-23.5 cm Circumference	(cA) 23.6-26.0 cm Circumference	(cA) 26.1-28.5 cm Circumference	(cA) 28.6-31.0 cm Circumference
Right <input type="checkbox"/> Beige Right <input type="checkbox"/> Black Left <input type="checkbox"/> Beige Left <input type="checkbox"/> Black	<input type="checkbox"/> 104416 <input type="checkbox"/> 104417 <input type="checkbox"/> 104406 <input type="checkbox"/> 104407	<input type="checkbox"/> 104414 <input type="checkbox"/> 104415 <input type="checkbox"/> 104404 <input type="checkbox"/> 104405	<input type="checkbox"/> 104412 <input type="checkbox"/> 104413 <input type="checkbox"/> 104402 <input type="checkbox"/> 104403	<input type="checkbox"/> 104418 <input type="checkbox"/> 104419 <input type="checkbox"/> 104408 <input type="checkbox"/> 104409	<input type="checkbox"/> 104420 <input type="checkbox"/> 104421 <input type="checkbox"/> 104410 <input type="checkbox"/> 104411
ReadyWrap® Extender Strap One size fits all					
<input type="checkbox"/> Beige <input type="checkbox"/> Black	<input type="checkbox"/> 104740 <input type="checkbox"/> 104741				
ReadyWrap® Liner Sold in pairs, available in black only					
	Small - Large		X-Large - XX-Large		
	Ankle Circumference (cB) 18-33 cm Calf Circumference (cC) 28-48 cm Groin Circumference (cG) 48-83 cm		Ankle Circumference (cB) 30-43 cm Calf Circumference (cC) 43-68 cm Groin Circumference (cG) 78-103 cm		
Below Knee <input type="checkbox"/> Black	<input type="checkbox"/> 149032		<input type="checkbox"/> 149033		
Thigh High <input type="checkbox"/> Black	<input type="checkbox"/> 149034		<input type="checkbox"/> 149035		