



People Health Care

ReadyWrap® Measuring below the knee only

For the full limb measuring form if oedema management above lower limb


Scan the QR code to access the form ►



Measuring form for Lower Limb - Complete one form per leg

Each garment comes with 1 x ReadyWrap® Liner

How to measure the leg: Measure the leg whilst patient is standing, if possible.



Write your measurement results in the box adjacent to measuring point

Please note

c = circumference, L = length

All length measurements should be longitudinal and should follow body contours.
All measurements to be in centimetres.

Patients Name: _____ Clinicians Name: _____


Date of measuring: _____ Contact Number: _____

Delivery: ☐ to orderer ☐ directly to patient




ReadyWrap® Toe

Circumference c / Length L	Small	Medium	Large	XL	XXL
cA around base of the toes	18-20.4cm	20.5-23.5cm	23.6-26cm	26.1-28.5cm	28.6-31cm



ReadyWrap® Foot CT

Circumference c / Length L	Small	Medium	Large	XL	XXL
cA around base of the toes	18-20.4cm	20.5-23.5cm	23.6-26cm	26.1-28.5cm	28.6-31cm
cA¹ widest part around the dorsum of the foot	20-25cm	24-29cm	28-33cm	32-37cm	36-41cm
cB 2cm above malleolus	18-23cm	20-28cm	25-33cm	30-38cm	35-43cm
LA-Y from the base of the little toe to the back of the heel	Up to 19cm – Regular			19-22cm – Long	



ReadyWrap® Calf

Circumference c / Length L	Small	Medium	Large	XL	XXL
cB 2cm above malleolus	18-20.4cm	20.5-23.5cm	23.6-26cm	26.1-28.5cm	28.6-31cm
cC widest part of calf	28-38cm	33-43cm	38-48cm	43-58cm	53-68cm
LD-B the length between point B (2cm above the malleolus) and 2cm below popliteal fossa	Up to 30cm – 30cm			Up to 35cm – 35cm	

Colour	Black	Beige
Left or Right Leg	Left	Right

ReadyWrap® Extender Strap

	Black	Beige
--	-------	-------

Additional ReadyWrap® Liners

Small – Large cB 18-33cm	XL – XXL cB 30-43cm
Sold in pairs, available in black only	

Correct ordering description example:

ReadyWrap® Foot CT Small Regular Black Left x1 / ReadyWrap® Calf Small 30cm (x1) Black

Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	XL	XXL
ReadyWrap® Toe					
Right ■ Beige	104416	104414	104412	104418	104420
Right ■ Black	104417	104415	104413	104419	104421
Left ■ Beige	104406	104404	104402	104408	104410
Left ■ Black	104407	104405	104403	104409	104411
ReadyWrap® Foot CT					
Right ■ Beige Reg up to 19cm	141727	141731	141735	141739	141743
Right ■ Black Reg up to 19cm	141729	141733	141737	141741	141745
Left ■ Beige Reg up to 19cm	141726	141730	141734	141738	141742
Left ■ Black Reg up to 19cm	141728	141732	141736	141740	141744
Right ■ Beige Long 19-22cm	141747	141751	141755	141759	141763
Right ■ Black Long 19-22cm	141749	141753	141757	141761	141765
Left ■ Beige Long 19-22cm	141746	141750	141754	141758	141762
Left ■ Black Long 19-22cm	141748	141752	141756	141760	141764
ReadyWrap® Calf					
■ Beige up to 30 cm	104450	104446	104442	104454	104784
■ Black up to 30 cm	104451	104447	104443	104455	104785
■ Beige up to 35 cm	104452	104448	104444	104456	104786
■ Black up to 35 cm	104453	104449	104445	104457	104787

ReadyWrap® Liner Sold in pairs, available in black only

	Small - Large	X-Large - XX-Large
Below Knee ■ Black	149032	149033
Thigh High ■ Black	149034	149035

ReadyWrap® Extender Strap One size fits all

■ Beige	104740
■ Black	104741

Scan to download
the **Compression
Selector App**



Scan here to access
**EPS and DAC
ordering support**



Pharmacy Support
Scan here to access
free pharmacy resources
and information



Scan here to access
our **FREE learning
modules and videos**



Scan here to contact
our **Customer
Services Team**

