

## ReadyWrap® Measuring the full lower limb only

For the below knee measuring form if oedema management below upper limb


Scan the QR code to access the form ▶



**Measuring form for full lower limb** - Complete one form per leg

Each garment comes with 1 x ReadyWrap® Liner

**How to measure the leg:** Measure the leg whilst patient is standing, if possible.



Write your measurement results in the box adjacent to measuring point

**Please note**

**c = circumference, L = length**

All length measurements should be longitudinal and should follow body contours.

All measurements to be in centimetres.

Patients Name: \_\_\_\_\_ Clinicians Name: \_\_\_\_\_

Date of measuring: \_\_\_\_\_ Contact Number: \_\_\_\_\_

Delivery: to orderer directly to patient

ReadyWrap® Thigh					
Circumference c / Length L	Small	Medium	Large	XL	XXL
cG 2cm below gluteal fold	48-63 cm	58-73 cm	68-83 cm	78-93 cm	88-103 cm
cE¹ 5cm above the knee	38-48cm	45-58cm	52-68cm	59-78cm	66-88cm
LG-E¹ the length between point cE¹ (5cm above the knee) and point cG (2cm below the gluteal fold)	Up to 25cm – 25cm			Up to 30cm – 30cm	

ReadyWrap® Knee					
Circumference c / Length L	Small	Medium	Large	XL	XXL
cD 2-3cm below popliteal fossa	28-38cm	33-43cm	38-48cm	43-58cm	53-68cm
cE¹ 5cm above the knee	38-48cm	45-58cm	52-68cm	59-78cm	66-88cm

ReadyWrap® Calf					
Circumference c / Length L	Small	Medium	Large	XL	XXL
cB 2cm above malleolus	18-23cm	20-28cm	25-33cm	30-38cm	35-43cm
cC widest part of calf	28-38cm	33-43cm	38-48cm	43-58cm	53-68cm
LD-B the length between point B (2cm above the malleolus) and 2cm below popliteal fossa	Up to 30cm – 30cm			Up to 35cm – 35cm	

ReadyWrap® Foot CT					
Circumference c / Length L	Small	Medium	Large	XL	XXL
cA around base of the toes	18-20.4cm	20.5-23.5cm	23.6-26cm	26.1-28.5cm	28.6-31cm
cA¹ widest part around the dorsum of the foot	20-25cm	24-29cm	28-33cm	32-37cm	36-41cm
cB 2cm above malleolus	18-23cm	20-28cm	25-33cm	30-38cm	35-43cm
LA-Y from the base of the little toe to the back of the heel	Up to 19cm – Regular			19-22cm – Long	

ReadyWrap® Toe					
Circumference c / Length L	Small	Medium	Large	XL	XXL
cA around base of the toes	18-20.4cm	20.5-23.5cm	23.6-26cm	26.1-28.5cm	28.6-31cm

<b>Colour</b>	Black	Beige
<b>Left or Right Leg</b>	Left	Right

<b>ReadyWrap® Extender Strap</b>	Black	Beige
<b>Additional ReadyWrap® Liners</b>	Small – Large cB 18-33cm	XL – XXL cB 30-43cm
Sold in pairs, available in black only		

*Correct ordering description example:*

*ReadyWrap® Foot CT Small Regular Black Left x1 / ReadyWrap® Calf Small 30cm (x1) Black*

Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	XL	XXL
<b>ReadyWrap® Toe</b>					
<b>Right</b> ■ Beige	104416	104414	104412	104418	104420
<b>Right</b> ■ Black	104417	104415	104413	104419	104421
<b>Left</b> ■ Beige	104406	104404	104402	104408	104410
<b>Left</b> ■ Black	104407	104405	104403	104409	104411
<b>ReadyWrap® Foot CT</b>					
<b>Right</b> ■ Beige Reg up to 19cm	141727	141731	141735	141739	141743
<b>Right</b> ■ Black Reg up to 19cm	141729	141733	141737	141741	141745
<b>Left</b> ■ Beige Reg up to 19cm	141726	141730	141734	141738	141742
<b>Left</b> ■ Black Reg up to 19cm	141728	141732	141736	141740	141744
<b>Right</b> ■ Beige Long 19-22cm	141747	141751	141755	141759	141763
<b>Right</b> ■ Black Long 19-22cm	141749	141753	141757	141761	141765
<b>Left</b> ■ Beige Long 19-22cm	141746	141750	141754	141758	141762
<b>Left</b> ■ Black Long 19-22cm	141748	141752	141756	141760	141764
<b>ReadyWrap® Calf</b>					
■ Beige up to 30 cm	104450	104446	104442	104454	104784
■ Black up to 30 cm	104451	104447	104443	104455	104785
■ Beige up to 35 cm	104452	104448	104444	104456	104786
■ Black up to 35 cm	104453	104449	104445	104457	104787
<b>ReadyWrap® Knee</b>					
■ Beige	104464	104462	104460	104466	104788
■ Black	104465	104463	104461	104467	104789
<b>ReadyWrap® Thigh</b>					
■ Beige up to 25 cm	104478	104474	104470	104482	104790
■ Black up to 25 cm	104479	104475	104471	104483	104791
■ Beige up to 30 cm	104480	104476	104472	104484	104792
■ Black up to 30 cm	104481	104477	104473	104485	104793

**ReadyWrap® Liner** Sold in pairs, available in black only

	Small - Large	X-Large - XX-Large
<b>Below Knee</b> ■ Black	149032	149033
<b>Thigh High</b> ■ Black	149034	149035

**ReadyWrap® Extender Strap** One size fits all

■ Beige	104740
■ Black	104741

Scan to download the **Compression Selector App**



Scan here to access **EPS and DAC** ordering support



**Pharmacy Support**  
Scan here to access free pharmacy resources and information



Scan here to access our **FREE learning modules and videos**



Scan here to contact our **Customer Services Team**

