

L&R's hand-finished **ActiLymph® Made-to-Measure** garments are manufactured using advanced yarn technology providing new levels of patient concordance. **ActiLymph®** has been developed specifically for lymphoedema and chronic oedema patients who want the quality and comfort of the popular **ActiLymph® European hosiery** in a custom fit garment.

Why choose ActiLymph® MTM?

- + A choice of 2 types of fabric to match your patients' needs
- + Available on FP10/GP10
- + Exceptional service for a hand-finished garment
- + Excellent patient feedback relating to fit and appearance
- + Class 1 (18-21mmHg)
Class 2 (23-32mmHg)
Class 3 (34-46mmHg)



"I have found the **ActiLymph® MTM Ease** garments incredible, the best I've worn!"
Feedback from a new ActiLymph® patient (August 2011)

* Delivery will be within 5 working days of receipt of a correctly completed order form; subject to circumstances beyond our control.
For further information on **ActiLymph® MTM** garments and training, please contact the L&R Customer Care Team on **08450 606707**
or visit our website at **www.Lohmann-Rauscher.co.uk**



- + Improved comfort and fit
- + Effective oedema containment
- + Easy application
- = Excellent patient concordance

Applied test norm and methods Norm/method
- ENV 12718 - NEN 8188 - RAL-GZ 387/1 - Öko-tex 100 - ISO 9001-2000



For **ActiLymph® MTM** delivery and distribution please contact Credenhill:
Email: sales@credenhill.co.uk
Fax: 0115 944 0437
Post: Credenhill Limited, 10 Cossall Industrial Estate, Ilkeston, Derbyshire, DE7 5UG.
Freephone: 0800 195 0650 or 0115 932 0144



Activa®, ActiLymph® and ActiGlide® are registered to L&R Medical UK Ltd.

M643 V2.2



ActiLymph®

MadeToMeasure
Flat Knit Hosiery



for Lymphoedema and Chronic Oedema

Available on Drug Tariff



5 working day service*



ActiLymph® MTM Ease



ActiLymph® MTM Facilitating patient comfort and concordance with a softer fabric - whilst delivering effective compression.

+ Available in the following classes and colours:

European Class			
	Class 1 (18-21mmHg)	Class 2 (23-32mmHg)	Class 3 (34-46mmHg)
Colour	Black ●	Black ●	Black ●
	Sand ●	Sand ●	Sand ●

+ Styles available:
 Below knee
 Thigh length
 Thigh length with waist attachment
 Two legged panty (includes cycle style pants)
 Open or closed toe

Your patients can now have extra comfort with ease...

- + Silky fabric that glides over the leg
- + Effortless application and easy removal
- + Discreet and elegant appearance
- + Ultra low profile seams
- + Fine toe seam and T-heel for ergonomic fit
- + Smooth grip top
- + Optional oblique toe



Application aids

Please note: With the MTM open toe stocking we include a pull-on sock to aid application. We also recommend the use of the **ActiGlide** application aid for all open and closed toe stockings. The **ActiGlide** is available on prescription. Contact your pharmacy for further information.



Pip Code: 284-5196

ActiLymph® MadeToMeasure

Style	Class	Identity Code
-------	-------	---------------

ActiLymph® MTM Ease

Below Knee A-D	Class 1	AL-EASEC1-AD
	Class 2	AL-EASEC2-AD
	Class 3	AL-EASEC3-AD

Thigh length A-G	Class 1	AL-EASEC1-AG
	Class 2	AL-EASEC2-AG
	Class 3	AL-EASEC3-AG

Thigh length with waist attachment A-GT	Class 1	AL-EASEC1-AGT
	Class 2	AL-EASEC2-AGT
	Class 3	AL-EASEC3-AGT

Two legged panty A-T (includes cycle style pants)	Class 1	AL-EASEC1-AT
	Class 2	AL-EASEC2-AT
	Class 3	AL-EASEC3-AT

ActiLymph® MTM Dura

Below Knee A-D	Class 2	AL-DURAC2-AD
	Class 3	AL-DURAC3-AD

Thigh length A-G	Class 2	AL-DURAC2-AG
	Class 3	AL-DURAC3-AG

Thigh length with waist attachment A-GT	Class 2	AL-DURAC2-AGT
	Class 3	AL-DURAC3-AGT

Two legged panty A-T (includes cycle style pants)	Class 2	AL-DURAC2-AT
	Class 3	AL-DURAC3-AT

Optional Extras

Style	Identity Code
Closed Toe	ALO - CT
2 Ankle pads (Dura only)	ALO - PAD
Silicone top band	ALO - BAND
Slip form	ALO - FORM
T-heel ankle seam	ALO - SEAM
Full Compression Body Part (not available in Class 1)	ALO - FCBODY
Oblique Toe (Slanted Toe)	ALO - OBT



ActiLymph® MTM Dura



ActiLymph® MTM Dura provides the fabric stiffness ideal for containing stubborn chronic oedema and lymphoedema, whilst maintaining patient comfort and concordance.

+ Available in the following classes and colours:

European Class			
	Class 1 (18-21mmHg)	Class 2 (23-32mmHg)	Class 3 (34-46mmHg)
Colour	Black ●	Black ●	Black ●
	Sand ●	Sand ●	Sand ●

+ Styles available:
 Below knee
 Thigh length
 Thigh length with waist attachment
 Two legged panty (includes cycle style pants)
 Open or closed toe

Performance and comfort in harmony...

- + For more stubborn and resistant chronic oedema
- + Low profile seams
- + Low profile toe seams and T-heel for ergonomic fit
- + Smooth grip top
- + Optional ankle pads
- + Optional oblique toe

