

ActiLymph®/Adore® European Class Hosiery - Patient Record

For patients with oedema This is not a made to measure form

Patient name _____

Date _____

Doctor/Clinician _____

Pharmacist _____

Clinic or Pharmacy Address:

(Please specify the address where the patient will collect their hosiery)

Request for ActiLymph®/Adore® European Class Hosiery

Hosiery measuring guide

Size	Small	Medium	Large	Extra large	Extra Extra large*
G Wide Topband	57 – 66cm	61 – 72cm	66 – 80cm	74 – 90cm	N/A
G Regular Topband/ Adore® Pattern Topband	48 – 56cm	50 – 60cm	53 – 66cm	57 – 73cm	61 – 80cm
C	31 – 39cm	33 – 41cm	35 – 44cm	39 – 48cm	44 – 52cm
B	19 – 21cm	21 – 23.5cm	23 – 27cm	27 – 32cm	32 – 37cm
LZ ActiLymph®	22 – 28cm	22 – 28cm	22 – 28cm	22 – 28cm	22 – 28cm
LZ Adore®	23 – 26cm	23 – 26cm	23 – 26cm	23 – 26cm	–

Style	Leg length (LD)	Foot length (LZ) ActiLymph®	Foot length (LZ) Adore®
Below knee	Petite <38cm	22 – 25cm	22 – 24cm
	Standard >38cm	22 – 28cm	23 – 26cm

* For availability of XXL garments, please see PIP codes on reverse.

Before applying hosiery, please read patient information leaflet enclosed within the box.

Garment prescribed

For guidance on which level of compression is most suitable for the patient's condition, please refer to the indications overleaf. **N.B.** Sizing charts vary from one manufacturer to another. Please ensure you prescribe products using the correct chart.

Size requested

Please state the number of packs (pairs) you require

Please note

If your patient has a long foot in proportion to their leg measurements, an open-toe option is advisable. Open toe is also appropriate where you require access to the toes.


If your patient's measurements do not fit into the stock sizes then a 'made to measure' garment may be more appropriate. For more information please see our website at www.Lohmann-Rauscher.co.uk or call customer care on 08450 606 707

Is the ActiGlide® Hosiery Applicator required?

(PIP Code: 284-5196)

Yes

No



Patient measurements

Hip	Left leg	Right leg
G: 5cm below groin		
C: Widest part of the calf		
B: Just above the malleolus (ankle bone)		
LZ: From back of the heel to the longest toe		
LD: Required for below knee garments only: Measure from the popliteal space to floor.		

Please tick a box to indicate garment required

	Colour	18-21mmHg Class 1	23-32mmHg Class 2	34-46mmHg Class 3	Hosiery Kit
ActiLymph® Thigh Length, Standard - Wide TopBand (Closed Toe)	Sand	<input type="checkbox"/>	<input type="checkbox"/>		
ActiLymph® Thigh Length, Standard - Wide TopBand (Open Toe)	Sand		<input type="checkbox"/>	<input type="checkbox"/>	
ActiLymph® Thigh Length, Standard - Regular TopBand (Closed Toe)	Black	<input type="checkbox"/>	<input type="checkbox"/>		
ActiLymph® Thigh Length, Standard - Regular TopBand (Open Toe)	Sand	<input type="checkbox"/>	<input type="checkbox"/>		
	Black	<input type="checkbox"/>	<input type="checkbox"/>		
Adore® Thigh Length, Standard - Patterned TopBand (Closed Toe)	Natural Tan	<input type="checkbox"/>	<input type="checkbox"/>		
Adore® Thigh Length, Standard - Patterned TopBand (Closed Toe)	Black	<input type="checkbox"/>	<input type="checkbox"/>		
Adore® Below Knee, Standard (Closed Toe)	Natural Tan	<input type="checkbox"/>	<input type="checkbox"/>		
Adore® Below Knee, Standard (Closed Toe)	Black	<input type="checkbox"/>	<input type="checkbox"/>		
Adore® Below Knee, Petite (Closed Toe)	Natural Tan	<input type="checkbox"/>	<input type="checkbox"/>		
Adore® Below Knee, Petite (Closed Toe)	Black	<input type="checkbox"/>	<input type="checkbox"/>		
ActiLymph® Below Knee, Standard (Closed Toe)	Sand	<input type="checkbox"/>	<input type="checkbox"/>		
	Black	<input type="checkbox"/>	<input type="checkbox"/>		
ActiLymph® Below Knee, Petite (Closed Toe)	Sand	<input type="checkbox"/>	<input type="checkbox"/>		
	Black	<input type="checkbox"/>	<input type="checkbox"/>		
ActiLymph® Below Knee, Standard (Open Toe)	Sand	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	Black	<input type="checkbox"/>	<input type="checkbox"/>		
ActiLymph® Below Knee, Petite (Open Toe)	Sand	<input type="checkbox"/>	<input type="checkbox"/>		
	Black	<input type="checkbox"/>	<input type="checkbox"/>		
ActiLymph® Hosiery Kit (1 below knee open toe stocking plus 2 compression liners)	Sand				<input type="checkbox"/>
	Black				<input type="checkbox"/>

ActiLymph® European Class Hosiery Range

(For Adore® range codes, please see Adore® code sheet)

Product & Style	Colour	Small	Medium	Large	Extra Large	Extra Extra Large
Class 1 (18-21mmHg)						
Provides light compression and should be used for early/mild chronic oedema or lymphoedema - where the oedema is light to moderate with little shape distortion. Suitable for chronic oedema, lymphoedema, lipoedema, prophylaxis, maintenance therapy and palliative use.						
Thigh Length Standard Wide TopBand	Sand	317-3614	317-3622	317-3630	317-3648	-
Thigh Length Standard Regular TopBand	Black	358-7474	358-7482	358-7490	358-7508	-
Thigh Length Open Toe Standard Regular TopBand	Sand	336-8842	336-8859	336-8867	336-8875	336-8883
Thigh Length Open Toe Standard Regular TopBand	Black	342-8257	342-8265	342-8273	342-8281	342-8299
Below Knee Standard	Sand	317-3531	317-3549	317-3556	317-3564	336-8834
Below Knee Standard	Black	342-8158	342-8166	342-8174	342-8182	342-8190
Below Knee Open Toe Standard	Sand	336-9402	336-9410	336-9428	336-9436	336-9444
Below Knee Open Toe Standard	Black	342-8208	342-8216	342-8224	342-8232	342-8240
Below Knee Petite	Sand	317-3572	317-3580	317-3598	317-3606	-
Class 2 (23-32mmHg)						
Provides medium compression and should be used in cases of chronic oedema or lymphoedema. ActiLymph® Class 2 hosiery can be used as part of the treatment for, and to provide maintenance therapy of chronic oedema and lymphoedema.						
Thigh Length Standard Wide TopBand	Sand	317-3739	317-3747	317-3754	317-3762	-
Thigh Length Open Toe Standard Wide TopBand	Sand	317-3853	317-3861	317-3879	317-3887	-
Thigh Length Standard Regular TopBand	Black	358-7516	358-7524	358-7532	358-7540	-
Thigh Length Open Toe Standard Regular TopBand	Sand	336-8990	336-9006	336-9014	336-9030	336-9048
Thigh Length Open Toe Standard Regular TopBand	Black	342-8398	342-8406	342-8414	342-8422	342-8430
Below Knee Standard	Sand	317-3655	317-3663	317-3671	317-3689	-
Below Knee Standard	Black	342-8307	342-8315	342-8323	342-8331	-
Below Knee Petite	Sand	317-3697	317-3705	317-3713	317-3721	-
Below Knee Open Toe Standard	Sand	317-3770	317-3788	317-3796	317-3804	336-9055
Below Knee Open Toe Standard	Black	342-8349	342-8356	342-8364	342-8372	342-8380
Below Knee Open Toe Petite	Sand	317-3812	317-3820	317-3838	317-3846	-
Class 3 (34-46mmHg)						
Provides strong compression and should be used for the maintenance of severe chronic oedema and lymphoedema, where resistant oedema features e.g. following a DVT or a history of reoccurring ulceration or where lymphatic damage is considerable. Use when lower compression has failed to control oedema.						
Thigh Length Open Toe Standard Wide TopBand	Sand	317-3978	317-3986	317-3994	317-4000	-
Below Knee Open Toe Standard	Sand	317-3895	317-3903	317-3911	317-3929	-
ActiLymph® Hosiery Kit						
The ActiLymph® Hosiery Kit is designed to be an alternative treatment to compression bandages for venous leg ulcer treatment. The hosiery kit facilitates self-care and enables carer involvement. The compression liner is applied first, thus allowing the Class 2 stocking to glide easily onto the limb. (Please note: 2 compression liners are included in the Kit - 1 to wear and 1 to be washed. Activa® 10mmHg liner packs are available separately, 3 per pack, in either white, sand or black.)						
Sand Kit (1 x Sand Class 2 Below Knee Open Toe Stocking + 2 x White Closed Toe Liners)		-	374-6914	374-6922	374-6930	374-6948
Black Kit (1 x Black Class 2 Below Knee Open Toe Stocking + 2 x Sand Closed Toe Liners)		-	387-0763	387-0771	387-0789	387-0797
ActiGlide® Hosiery Applicator						
ActiGlide® Hosiery Applicator	One Size Fits All	284-5196	-	-	-	-

For further information please call **08450 606 707**
to speak directly with a Customer Care Advisor
or visit our website www.Lohmann-Rauscher.co.uk

