



People.Health.Care.

# ActiLymph<sup>®</sup> Made To Measure

## Flat knit hosiery for Lymphoedema and Chronic Oedema



### ActiLymph<sup>®</sup> MTM Ease

Available in European Class 1 (18-21mmHg), Class 2 (23-32mmHg) and Class 3 (34-46mmHg).

**ActiLymph<sup>®</sup> MTM Ease** is silky, effective, comfortable to wear and is ideal for most forms of chronic oedema and lymphoedema.



### ActiLymph<sup>®</sup> MTM Dura

Available in European Class 2 (23-32mmHg) and Class 3 (34-46mmHg).

**ActiLymph<sup>®</sup> MTM Dura** is designed for more difficult to manage and resistant chronic oedema and lymphoedema.

5 working day service\*

# ActiLymph<sup>®</sup> Made To Measure

## Flat knit hosiery for Lymphoedema and Chronic Oedema

### How to measure the leg:

Ask the patient to stand, if possible.

Mark and measure the leg as follows:-

- cG 2cm below Gluteal Fold for circumference
- IG 2cm below Gluteal Fold to ground for length
- cF Mid Thigh for circumference
- IF Mid Thigh to ground for length
- cE Across middle of Kneecap for circumference
- IE Middle of Kneecap to ground for length
- cD 2cm below Popliteal Fossa for circumference
- ID 2cm below Popliteal Fossa to ground for length
- cC Widest part of Calf for circumference
- IC Widest part of Calf to ground for length
- cB1 Where Calf starts for circumference
- IB1 Where Calf starts to ground for length
- cB 2cm Above Malleolus for circumference
- IB 2cm Above Malleolus to ground for length
- cY Around Instep and Heel with foot dorsi-flexed for circumference
- cA Around base of Toes for circumference
- IA Heel to base of Toes for length
- IZ Heel to end of Large Toe (Total length of foot)
- IA1 From the Heel to the base of Large Toe
- IA2 From the Heel to the base of Small Toe

**NB** c = circumference, l = length

Unless specified, all length measurements should be longitudinal and should not follow body contours

### How to measure the panty:

Ask patient to put their hands on their waistline:-

Measure (cT) waist circumference

Measure (IT) waist to ground

Measure (cH) greatest circumference around hips

Measure (IH) length from hip to ground

### Measure the following lengths along contours of body:

Measure (IGT) length from waist to 2cm below gluteal fold

Measure (IK1T) length in the centre at the front, from waist to crotch while applying gentle pressure and following contours

Measure (IK2T) length at the back from the waist over the buttocks to the gluteal fold while applying gentle pressure and following contours

\*Delivery will be within 5 working days of receipt of a correctly completed order form; subject to circumstances beyond our control.

**Please send the completed Order Form to Credenhill via:**

Email: [sales@credenhill.co.uk](mailto:sales@credenhill.co.uk)

Fax: **0115 944 0437**

Post: **Credenhill Limited, 10 Cossall Industrial Estate, Ilkeston, Derbyshire, DE7 5UG.**

Freephone: **0800 195 0650** or **0115 932 0144**

# ActiLymph® Made To Measure Order Form

**1** Name of Pharmacy or Clinic: \_\_\_\_\_  
 Contact name: \_\_\_\_\_  
 Account Number (if known): \_\_\_\_\_  
 Branch Number (if applicable): \_\_\_\_\_  
 Address: \_\_\_\_\_  
 \_\_\_\_\_  
 Postcode: \_\_\_\_\_  
 Tel: \_\_\_\_\_  
 Practitioner name: \_\_\_\_\_ Tel: \_\_\_\_\_  
 Patient name: \_\_\_\_\_  
 Order Ref. or Prescription No.: \_\_\_\_\_

Please complete this form and send to Credenhill Limited.  
 Fax: 0115 944 0437, or email: sales@credenhill.co.uk

**2** **Left leg** **Right leg**  
**Length Circumference Circumference Length**

**IT cT cT IT**  
**IH cH cH IH**  
**IG cG cG IG**  
**IF cF cF IF**  
**IE cE cE IE**  
**ID cD cD ID**  
**IC cC cC IC**  
**IB<sup>1</sup> cB<sup>1</sup> cB<sup>1</sup> IB<sup>1</sup>**  
**IB cB cB IB**  
**cY cY**  
**cA cA**  
**IA IA**  
**IZ IZ**

**NB** All length measurements should be longitudinal and should not follow body contours. LENGTH MEASUREMENTS TAKEN FROM POINTS ON THE LEG SHOULD BE MEASURED VERTICALLY TO THE GROUND. All measurements to be in centimetres.

**3** **Panty**  
**NB** Measure all lengths along contours of body. All measurements to be in centimetres.

**4** **Ankle Pads**  
**NB** All length measurements should be longitudinal and should not follow foot contours. All measurements to be in centimetres.

**LMA:** from the centre of the medial malleolus, in horizontal line over the Achilles tendon to the middle of the lateral malleolus.  
**LMO:** from the centre of the medial malleolus, in vertical line under the sole of the foot to the middle of the lateral malleolus.

**\* Left Foot** **\* Right Foot**  
**LMA** \_\_\_\_\_ **LMA** \_\_\_\_\_  
**LMO** \_\_\_\_\_ **LMO** \_\_\_\_\_

If Ankle pads are requested please select which you require (available in Dura only):  
**Left**  Medial  Lateral  
**Right**  Medial  Lateral

**5** **Oblique Toe (Slanted Toe)**  
**NB** All length measurements should be longitudinal and should not follow foot contours. All measurements to be in centimetres.

**Left Foot** **Right Foot**  
**IA<sup>1</sup>** **IA<sup>1</sup>**  
**IA<sup>2</sup>** **IA<sup>2</sup>**

## 6 Main Garment Details

Quantity required	Garment required		Colour required
	<b>ActiLymph® MTM Ease</b> <input type="checkbox"/> Class 1 <input type="checkbox"/> Class 2 <input type="checkbox"/> Class 3		<b>ActiLymph® MTM Dura</b> <input type="checkbox"/> Class 2 <input type="checkbox"/> Class 3
	<input type="checkbox"/> Black <input type="checkbox"/> Sand		
<b>Style</b>			
<input type="checkbox"/> Below Knee (A-D)	<input type="checkbox"/> Thigh Length (A-G)	<input type="checkbox"/> Thigh Length with Waist Attachment (A-GT)	<input type="checkbox"/> 2 Legged Panty Hose (A-T)
<b>Options</b>			
<input type="checkbox"/> Closed Toe (measurements IA & IZ required) <input type="checkbox"/> Open Toe (IA only) <input type="checkbox"/> Closed Oblique (Slanted) Toe** (measurement cA, IA1, IA2 & IZ required)			
<input type="checkbox"/> T Heel Ankle Seam <input type="checkbox"/> Open Oblique (Slanted) Toe** (measurements cA, IA1 & IA2 required) <input type="checkbox"/> Ankle pads* (measurements LMO & LMA required- indicate requirements in section 4)			
<input type="checkbox"/> 3cm Silicone Top Band (only below knee) <input type="checkbox"/> 5cm Silicone Top Band <input type="checkbox"/> Slip Form <input type="checkbox"/> Full Compression Body Part (Class 2&3 only. measurements IK1T, IK2T & IGT required)			

**7** Product Code (see Drug Tariff Product Code Guide on reverse of form for clarification on code required): \_\_\_\_\_  
 \_\_\_\_\_  
 Additional Comments: \_\_\_\_\_

**NB** Submission of this form confirms the patient's authorisation for the information to be retained for the production and supply of their garment. L&R or Credenhill cannot be held responsible for incorrect measurements. To avoid delay, please ensure all appropriate measurements are completed and checked.

# Actilymph® Made To Measure Drug Tariff Codes



in partnership with



**Your prescription must include:**

- Step 1** Style/Class code of main garment
- Step 2** The codes for each optional extra required
- Step 3** The quantity of garments required
- Step 4** The colour

e.g. for 1 Sand below knee Ease stocking in class 3, with closed toes and T Heel ankle seam, the prescription must include: 1x **AL-EASEC3-AD** together with optional extras 1x **ALO-CT** and 1x **ALO-SEAM, Sand**.

Please tick where applicable to indicate **Main Garment** choice and each **Optional Extra** included. Items listed below have Drug Tariff coding and associated prices. The Drug Tariff codes are required on prescriptions for **Actilymph® MTM**.

Main Garments				Optional Extras						
Style	Class	Main Garment	Colours	Closed Toe	Oblique Toe	2 Ankle Pads <small>Indicate Medial or Lateral option overleaf</small>	Silicone Topband <small>Indicate 3cm or 5cm option overleaf</small>	Slip Form	T Heel Ankle Seam	Full Compression Body Part
<input type="checkbox"/> <b>Actilymph® MTM Ease</b> Below Knee A-D	<input type="checkbox"/> Class 1	<input type="checkbox"/> AL-EASEC1-AD	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT		<input type="checkbox"/> ALO-BAND		<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 2	<input type="checkbox"/> AL-EASEC2-AD	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT		<input type="checkbox"/> ALO-BAND		<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 3	<input type="checkbox"/> AL-EASEC3-AD	<input type="checkbox"/> Sand <input type="checkbox"/> Black*	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT		<input type="checkbox"/> ALO-BAND		<input type="checkbox"/> ALO-SEAM	
<input type="checkbox"/> <b>Actilymph® MTM Ease</b> Thigh Length A-G	<input type="checkbox"/> Class 1	<input type="checkbox"/> AL-EASEC1-AG	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT		<input type="checkbox"/> ALO-BAND	<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 2	<input type="checkbox"/> AL-EASEC2-AG	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT		<input type="checkbox"/> ALO-BAND	<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 3	<input type="checkbox"/> AL-EASEC3-AG	<input type="checkbox"/> Sand <input type="checkbox"/> Black*	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT		<input type="checkbox"/> ALO-BAND	<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
<input type="checkbox"/> <b>Actilymph® MTM Ease</b> Thigh Length with waist attachment A-GT	<input type="checkbox"/> Class 1	<input type="checkbox"/> AL-EASEC1-AGT	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT			<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 2	<input type="checkbox"/> AL-EASEC2-AGT	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT			<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 3	<input type="checkbox"/> AL-EASEC3-AGT	<input type="checkbox"/> Sand <input type="checkbox"/> Black*	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT			<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
<input type="checkbox"/> <b>Actilymph® MTM Ease</b> Two legged panty A-T (includes cycle style pants)	<input type="checkbox"/> Class 1	<input type="checkbox"/> AL-EASEC1-AT	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT				<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 2	<input type="checkbox"/> AL-EASEC2-AT	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT				<input type="checkbox"/> ALO-SEAM	<input type="checkbox"/> ALO-FCBODY
	<input type="checkbox"/> Class 3	<input type="checkbox"/> AL-EASEC3-AT	<input type="checkbox"/> Sand <input type="checkbox"/> Black*	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT				<input type="checkbox"/> ALO-SEAM	<input type="checkbox"/> ALO-FCBODY
<input type="checkbox"/> <b>Actilymph® MTM Dura</b> Below Knee A-D	<input type="checkbox"/> Class 2	<input type="checkbox"/> AL-DURAC2-AD	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-PAD	<input type="checkbox"/> ALO-BAND		<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 3	<input type="checkbox"/> AL-DURAC3-AD	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-PAD	<input type="checkbox"/> ALO-BAND		<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 2	<input type="checkbox"/> AL-DURAC2-AG	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-PAD	<input type="checkbox"/> ALO-BAND	<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
<input type="checkbox"/> <b>Actilymph® MTM Dura</b> Thigh Length A-G	<input type="checkbox"/> Class 3	<input type="checkbox"/> AL-DURAC3-AG	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-PAD	<input type="checkbox"/> ALO-BAND	<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 2	<input type="checkbox"/> AL-DURAC2-AGT	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-PAD		<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
	<input type="checkbox"/> Class 3	<input type="checkbox"/> AL-DURAC3-AGT	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-PAD		<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	
<input type="checkbox"/> <b>Actilymph® MTM Ease</b> Two legged panty A-T (includes cycle style pants)	<input type="checkbox"/> Class 2	<input type="checkbox"/> AL-DURAC2-AT	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-PAD		<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	<input type="checkbox"/> ALO-FCBODY
	<input type="checkbox"/> Class 3	<input type="checkbox"/> AL-DURAC3-AT	<input type="checkbox"/> Sand <input type="checkbox"/> Black	<input type="checkbox"/> ALO-CT <input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-0BT	<input type="checkbox"/> ALO-PAD		<input type="checkbox"/> ALO-FORM	<input type="checkbox"/> ALO-SEAM	<input type="checkbox"/> ALO-FCBODY

\* Delivery will be within 5 working days of receipt of a correctly completed order form; subject to circumstances beyond our control.

