



People.Health.Care.

## Activa® British Standard Hosiery - Patient Record

For patients **without oedema** This is not a made to measure form

Patient name \_\_\_\_\_

Date \_\_\_\_\_

Doctor/Clinician \_\_\_\_\_

Pharmacist \_\_\_\_\_

### Clinic or Pharmacy Address:

(Please specify the address where the patient will collect their hosiery)

## Request for Activa® British Standard hosiery on FP10/GP10

### Hosiery measuring guide

A size indicator can be found by a coloured thread in the stocking

	Red	Blue	Green	Yellow	Purple
Size	Small	Medium	Large	Extra large	Extra Extra large*
Hip	Up to 107cm	107 – 112cm	112 – 122cm	122 – 132cm	N/A
A	40.5 – 53.0cm	43.0 – 55.0cm	48.0 – 60.5cm	53.0 – 66.0cm	N/A
B	30.5 – 37.5cm	33.0 – 40.0 cm	35.5 – 43.0cm	38.0 – 46.0cm	40.0 – 49.0cm
C	19.5 – 25.5cm	21.5 – 27.5cm	22.5 – 29.5cm	23.0 – 32.0cm	23.5 – 34.0cm
D	20.5 – 24.0cm	23.0 – 26.0cm	25.5 – 29.5cm	26.5 – 32.5cm	26.5 – 32.5cm


\* For availability of XXL garments, please see PIP codes on reverse.

### Garment prescribed

For guidance on which level of compression is most suitable for the patient's condition, please refer to the indications overleaf. **N.B.** Sizing charts vary from one manufacturer to another. Please ensure you prescribe products using the correct chart.

Size requested

Please state the number of packs you require



### Patient measurements

HIP	Left leg	Right leg
Hip		
A: 5cm below groin		
B: Widest part of the calf		
C: Just above the malleolus (ankle bone)		
D: From back of heel to the longest toe		

### Please note

If your patient has a long foot in proportion to their other leg measurements, an open-toe option is advisable. Open toe is also appropriate where you require access to the toes.

If your patient's measurements do not fit into the stock sizes then a 'made to measure' garment may be more appropriate. If oedema is present in the patient's limb, you may choose to consider applying ActiLymph® European Compression Hosiery. For more information please see our website at [www.Lohmann-Rauscher.co.uk](http://www.Lohmann-Rauscher.co.uk) or call customer care on 08450 606 707

### Is the ActiGlide® Hosiery Applicator required?

(PIP Code: 284-5196) Yes  No

Please tick a box to indicate garment required	Colour	10mmHg	14-17mmHg Class 1	18-24mmHg Class 2	25-35mmHg Class 3	40mmHg
Compression Liner (Closed Toe)†	White	<input type="checkbox"/>				
	Sand	<input type="checkbox"/>				
	Black	<input type="checkbox"/>				
Compression Liner (Open Toe)	Sand	<input type="checkbox"/>				
Below Knee Stockings (Closed Toe)	Sand		<input type="checkbox"/>	<input type="checkbox"/>		
	Honey		<input type="checkbox"/>	<input type="checkbox"/>		
	Black		<input type="checkbox"/>	<input type="checkbox"/>		
Below Knee Stockings (Open Toe)	Sand		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Thigh Length Stockings (Closed Toe)	Sand		<input type="checkbox"/>	<input type="checkbox"/>		
	Honey		<input type="checkbox"/>	<input type="checkbox"/>		
	Black		<input type="checkbox"/>	<input type="checkbox"/>		
Thigh Length Stockings (Open Toe)	Sand			<input type="checkbox"/>	<input type="checkbox"/>	
Unisex Sock	Black		<input type="checkbox"/>	<input type="checkbox"/>		
	Brown			<input type="checkbox"/>		
Unisex Patterned Sock	Black		<input type="checkbox"/>	<input type="checkbox"/>		
	Black		<input type="checkbox"/>			
Tights (Tights are not available on prescription)	Natural			<input type="checkbox"/>		
Anklet (See website for sizing chart)	Sand			<input type="checkbox"/>		
Kneecap (See website for sizing chart)	Sand			<input type="checkbox"/>		
Leg Ulcer Hosiery Kit (1 x below knee open toe stocking plus 2 compression liners)	Sand Stocking/White Liners					<input type="checkbox"/>
	Black Stocking/White Liners					<input type="checkbox"/>

† Search under 'Activa Leg Ulcer Hosiery Kit' to find Activa Compression Liners on your ordering system.

# Activa® British Standard Hosiery Range

Product & Style Colour of size indicator thread in stocking	Colour	Small Red	Medium Blue	Large Green	Extra Large Yellow	Extra Extra Large Purple
<b>Liner Pack - 10mmHg</b>						
The Activa® liner provides light compression for maintenance therapy, palliative care and for patients who are non-concordant with higher levels of compression, working on the assumption that some compression is better than none. They may be used with other compression stockings to increase compression on the limb, or to build the compression level in layers if the patient struggles to apply a higher class of hosiery.						
Stocking Liner (3 Singles)	White	311-8767	311-8775	311-8783	311-8791	311-8809
Stocking Liner (3 Singles)	Sand	326-9750	326-9768	326-9776	326-9784	326-9792
Stocking Liner (3 Singles)	Black	387-0805	387-0813	387-0821	387-0839	387-0847
Stocking Liner Open Toe (3 Singles)	Sand	340-3557	340-3565	340-3573	340-3581	340-3599
<b>Class 1 - 14-17mmHg</b>						
Activa® Class 1, medically proven, graduated support hosiery provides effective compression in the treatment of superficial and early varicose veins, including those during pregnancy. Class 1 hosiery can reduce the risk of Deep Vein Thrombosis in lower risk groups during long-distance travel or in long periods of inactivity.						
Thigh Length	Sand	259-0305	259-0339	259-0362	259-0396	-
Thigh Length	Honey	259-0313	259-0347	259-0370	259-0404	-
Thigh Length	Black	259-0321	259-0354	259-0388	259-0412	-
Below Knee	Sand	259-0420	259-0446	259-0461	259-0487	-
Below Knee	Honey	259-0438	259-0453	259-0479	259-0495	-
Below Knee	Black	278-2373	278-2381	278-2399	278-2407	-
Below Knee Open Toe	Sand	268-1500	268-1518	268-1526	268-1534	317-2723
Unisex Sock	Black	278-2415	278-2423	278-2431	278-2449	-
Unisex Patterned Sock	Black	356-9993	357-0009	357-0017	357-0025	-
Tights	Black	286-9998	287-0004	287-0012	287-0020	-
<b>Class 2 - 18-24mmHg</b>						
Activa® Class 2, medically approved, graduated support hosiery provides compression for medium severity varicose veins or prevention of venous leg ulcers. Class 2 hosiery can reduce the risk of Deep Vein Thrombosis in higher risk groups during long distance travel.						
Thigh Length	Sand	259-0503	259-0537	259-0560	259-0594	-
Thigh Length	Honey	259-0511	259-0545	259-0578	259-0602	-
Thigh Length	Black	259-0529	259-0552	259-0586	259-0610	-
Thigh Length Open Toe	Sand	259-0859	259-0867	259-0875	259-0883	-
Below Knee	Sand	259-0628	259-0685	259-0743	259-0784	-
Below Knee	Honey	259-0651	259-0701	259-0750	259-0792	-
Below Knee	Black	259-0669	259-0735	259-0776	259-0800	-
Below Knee Open Toe	Sand	259-0818	259-0826	259-0834	259-0842	317-2707
Unisex Sock	Black	259-4463	259-0891	259-0917	259-0933	-
Unisex Sock	Brown	259-4471	259-0909	259-0925	259-0941	-
Unisex Patterned Sock	Black	357-0033	357-0041	357-0058	357-0066	-
Anklets	Sand	259-0958	259-0966	259-0974	259-0982	-
Kneecaps	Sand	259-1204	259-1212	259-1220	259-1238	-
Tights	Natural	287-0038	287-0046	287-0053	287-0061	-
<b>Class 3 - 25-35mmHg</b>						
Activa® Class 3 hosiery provides effective management of chronic venous insufficiency and eases away the pain and discomfort of severe varicose veins and chronic leg ulcers. The hosiery can also be used in the prevention of leg ulcer recurrence.						
Thigh Length Open Toe	Sand	259-1287	259-1295	259-1303	259-1311	-
Below Knee Open Toe	Sand	259-1246	259-1253	259-1261	259-1279	310-2977
<b>40mmHg Leg Ulcer Hosiery Kit</b>						
The Activa® Leg Ulcer Hosiery Kit is designed to be an alternative treatment and maintenance solution in the management of uncomplicated venous leg disorders. The hosiery kit is suitable for patients who will not accept compression bandages and is particularly suitable for self-caring and younger patients.						
Sand (1 Sand Class 3 Stocking + 2 White Liners)	-	311-9021	311-9047	311-9054	311-9070	311-9088
Black (1 Black Class 3 Stocking + 2 Sand Liners)	-	366-5593	366-5601	366-5619	366-5627	366-5635
<b>ActiGlide® Hosiery Applicator</b>						
ActiGlide® Hosiery Applicator	One Size Fits All	284-5196				



For further information please call **08450 606 707**  
to speak directly with a Customer Care Advisor  
or visit our website [www.Lohmann-Rauscher.co.uk](http://www.Lohmann-Rauscher.co.uk)

Activa®, ActiLymph® and ActiGlide® are registered to L&R Medical UK Ltd, formerly known as Activa Healthcare Ltd



M002 V4.3