


# ReadyWrap®

## Measuring for below knee garments.

**How to measure the leg:** Ask the patient to stand, if possible. Mark and measure the leg as follows



Write your measurement results in the **box** adjacent to measuring point

### Circumference points:

- cC** widest part of Calf
- cB** 2 cm Above Malleolus
- cY** around Instep and Heel with foot dorsi-flexed
- cA** around base of Toes

### Lengths:

- IA-Y** heel to base of Toes
- ID-B** the length between point B (2 cm above the malleolus) and point D (2 cm below popliteal fossa)

### NB

**c** = circumference, **I** = length  
 All length measurements should be longitudinal and should follow body contours.  
 All measurements to be in centimetres.

Delivery:  to orderer  directly to patient

Name and address of recipient:

Patient Reference \_\_\_\_\_

Name/Last Name \_\_\_\_\_

Clinicians Name \_\_\_\_\_

Street \_\_\_\_\_

Clinicians Work Base \_\_\_\_\_

Postal Code/Place \_\_\_\_\_

Phone \_\_\_\_\_

Phone \_\_\_\_\_

Date of Measurement \_\_\_\_\_

Order Date \_\_\_\_\_

# Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	Extra Large	Extra Extra Large
<b>ReadyWrap® Calf</b>					
(I-D-B) (2 x lengths)	Calf Circumference (cC) 28-38 cm Ankle Circumference (cE) 18-23 cm	Calf Circumference (cC) 33-43 cm Ankle Circumference (cE) 20-28 cm	Calf Circumference (cC) 38-48 cm Ankle Circumference (cE) 25-33 cm	Calf Circumference (cC) 43-58 cm Ankle Circumference (cE) 30-38 cm	Calf Circumference (cC) 53-68 cm Ankle Circumference (cE) 35-43 cm
■ Beige up to 30 cm ■ Black up to 30 cm ■ Beige up to 35 cm ■ Black up to 35 cm	104450 104451 104452 104453	104446 104447 104448 104449	104442 104443 104444 104445	104454 104455 104456 104457	104784 104785 104786 104787
<b>ReadyWrap® Foot</b>					
(Side)	Foot Length (I-A-Y) up to 21 cm Ankle Circumference (cE) 18-23 cm Circumference (cY) 28-33 cm	Foot Length (I-A-Y) up to 23 cm Ankle Circumference (cE) 20-28 cm Circumference (cY) 30-38 cm	Foot Length (I-A-Y) up to 26 cm Ankle Circumference (cE) 25-33 cm Circumference (cY) 35-43 cm	Foot Length (I-A-Y) up to 29 cm Ankle Circumference (cE) 30-38 cm Circumference (cY) 40-48 cm	
Right ■ Beige Right ■ Black Left ■ Beige Left ■ Black	104436 104437 104428 104429	104434 104435 104426 104427	104432 104433 104424 104425	104438 104439 104430 104431	
<b>ReadyWrap® Toes</b>					
(cA)	(cA) 18.0-20.4 cm Circumference	(cA) 20.5-23.5 cm Circumference	(cA) 23.6-26.0 cm Circumference	(cA) 26.1-28.5 cm Circumference	(cA) 28.6-31.0 cm Circumference
Right ■ Beige Right ■ Black Left ■ Beige Left ■ Black	104416 104417 104406 104407	104414 104415 104404 104405	104412 104413 104402 104403	104418 104419 104408 104409	104420 104421 104410 104411
<b>ReadyWrap® Extender Strap</b>					
	One size fits all				
■ Beige ■ Black	104740 104741				

Please mark (tick or cross) where applicable to indicate the garment choice and then enter the sizing, quantity and reference code.

Garment (Left Leg)	Length/Side	Colour	Size	Quantity	Code	Garment (Right Leg)	Length/Side	Colour	Size	Quantity	Code
Calf	30 cm 35 cm	Beige Black				Calf	30 cm 35 cm	Beige Black			
Foot	Left	Beige Black				Foot	Right	Beige Black			
Toe	Left	Beige Black				Toe	Right	Beige Black			
Extender Strap		Beige Black				Extender Strap		Beige Black			

## Customer service

Lohmann & Rauscher offers you an easy solution for self-care compression. If you have any questions or require additional information, please contact our customer solutions team, from Monday to Friday:

Tel.: 08450 606 707 / 44 (0)1283 576800

Fax: 01283 576808

Email: customersolutions@uk.lmed.com

