



People.Health.Care.

ReadyWrap™

Arm sizing chart and order form.



Now on Drug Tariff
XL arm sleeve size





People.Health.Care.

ReadyWrap™

How to measure the arm:

Measure the arm in a supported and relaxed position with the elbow slightly flexed. The measurement points refer to the following anatomical landmarks:

Circumference points:

- cA Widest width of the palm, just below the metacarpal phalangeal joints
- cC Narrowest point of the wrist
- cE Around the elbow joint and into the crease of the elbow
- cG Top of the arm, just below the axillary crease

Lengths:

- IA-C the length between point A (widest width of the palm, just below the metacarpal phalangeal joints) and point C (narrowest point of the wrist)
- IC-G the length between point C (narrowest point of the wrist) and point G (top of the arm, just below the axillary crease)

NB c = circumference, l = length

Measurements required for each ReadyWrap upper limb garments:

For the Armsleeve Garment:

cG, cE, cC, IC-G

For the Gauntlet Garment:

cC, cA, IA-C

For further information on ReadyWrap garments and training please contact L&R Customer Care on **08450 606707** and visit our website **www.Lohmann-Rauscher.co.uk/readywrap** or contact your local L&R representative.

Delivery will be within 7 working days of receipt of a correctly completed order form; subject to circumstances beyond our control.

Please send the completed order form to Credenhill via:

Email: **sales@credenhill.co.uk**

Fax: **0115 944 0437**

Post: **Credenhill Limited, 10 Cossall Industrial Estate, Ilkeston, Derbyshire, DE7 5UG.**

Phone: **0115 932 0144**





People.Health.Care.

ReadyWrap™ Order Form

Name of Pharmacy or Clinic: _____

Contact Name: _____

Account Number (if known): _____ Branch Number (if applicable): _____

Address: _____

Postcode: _____

Tel: _____ Email: _____


Practitioner Name: _____ Tel: _____

Order Ref. or Prescription No.: _____

Patient Name/Ref: _____

Date of Measurement: ____/____/____ Order Date: ____/____/____

NB All length measurements should follow body contour. All measurements to be in centimetres.



Garment (Armsleeve)	Length/Side	Colour	Size	Quantity	Code <small>transfer code from back of sheet</small>
Arm	43cm <input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/>			
	46cm <input type="checkbox"/>				
	49cm <input type="checkbox"/>				
Gauntlet	Right <input type="checkbox"/> Left <input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/>			



People.Health.Care.

ReadyWrap™ Order Form

Please tick where applicable to indicate your garment choice

ReadyWrap™ Armsleeve (IC-G) (3 x lengths)	Small Axilla Circumference (cG) 22-29cm Elbow Circumference (cE) 19-26cm Wrist Circumference (cC) 13.5-17cm	Medium Axilla Circumference (cG) 26-36cm Elbow Circumference (cE) 24-33cm Wrist Circumference (cC) 15.5-20cm	Large Axilla Circumference (cG) 31-43cm Elbow Circumference (cE) 29-40cm Wrist Circumference (cC) 17.5-23cm	Extra Large Axilla Circumference (cG) 36-50cm Elbow Circumference (cE) 30-43cm Wrist Circumference (cC) 17.5-23cm
Short <input type="checkbox"/> Beige up to 43cm <input type="checkbox"/> Black up to 43cm Average <input type="checkbox"/> Beige up to 46cm <input type="checkbox"/> Black up to 46cm Tall <input type="checkbox"/> Beige up to 49.5cm <input type="checkbox"/> Black up to 49.5cm	<input type="checkbox"/> 104514 <input type="checkbox"/> 104515 <input type="checkbox"/> 104516 <input type="checkbox"/> 104517 <input type="checkbox"/> 104518 <input type="checkbox"/> 104519	<input type="checkbox"/> 104508 <input type="checkbox"/> 104509 <input type="checkbox"/> 104510 <input type="checkbox"/> 104511 <input type="checkbox"/> 104512 <input type="checkbox"/> 104513	<input type="checkbox"/> 104502 <input type="checkbox"/> 104503 <input type="checkbox"/> 104504 <input type="checkbox"/> 104505 <input type="checkbox"/> 104506 <input type="checkbox"/> 104507	<input type="checkbox"/> 104520 <input type="checkbox"/> 104521 <input type="checkbox"/> 104522 <input type="checkbox"/> 104523 <input type="checkbox"/> 104524 <input type="checkbox"/> 104525
ReadyWrap™ Gauntlet	Small Length (IA-C) up to 12.5cm Wrist Circumference (cC) 13.5-17cm Palm Circumference (cA) 18-20cm	Medium Length (IA-C) up to 13.5cm Wrist Circumference (cC) 15.5-20cm Palm Circumference (cA) 20-22cm	Large Length (IA-C) up to 14.5cm Wrist Circumference (cC) 17.5-23cm Palm Circumference (cA) 22-24cm	Extra Large Length (IA-C) up to 15.5cm Wrist Circumference (cC) 19.5-26cm Palm Circumference (cA) 24-26cm
Right <input type="checkbox"/> Beige Right <input type="checkbox"/> Black Left <input type="checkbox"/> Beige Left <input type="checkbox"/> Black	<input type="checkbox"/> 104498 <input type="checkbox"/> 104499 <input type="checkbox"/> 104490 <input type="checkbox"/> 104491	<input type="checkbox"/> 104496 <input type="checkbox"/> 104497 <input type="checkbox"/> 104488 <input type="checkbox"/> 104489	<input type="checkbox"/> 104494 <input type="checkbox"/> 104495 <input type="checkbox"/> 104486 <input type="checkbox"/> 104487	<input type="checkbox"/> 104500 <input type="checkbox"/> 104501 <input type="checkbox"/> 104492 <input type="checkbox"/> 104493

Delivery will be within 7 working days of receipt of a correctly completed order form; subject to circumstances beyond our control.

Please send the completed order form to Credenhill via:

Email: sales@credenhill.co.uk

Fax: **0115 944 0437**

Post: **Credenhill Limited, 10 Cossall Industrial Estate, Ilkeston, Derbyshire, DE7 5UG.**

Phone: **0115 932 0144**