


# ReadyWrap®

## Measuring the full lower limb.

**How to measure the leg:** Ask the patient to stand, if possible. Mark and measure the leg as follows



Write your measurement results in the **box** adjacent to measuring point

### Circumference points:

- cG** 2 cm below Gluteal Fold
- cE'** 5 cm above the Knee
- cD** 2-3 cm below Popliteal Fossa
- cC** widest part of Calf
- cB** 2 cm Above Malleolus
- cY** around Instep and Heel with foot dorsi-flexed
- cA** around base of Toes

### Lengths:

- IA-Y** heel to base of Toes
- ID-B** the length between point B (2 cm above the malleolus) and point D (2 cm below popliteal fossa)
- IG-E'** the length between point E' (5 cm above the knee) and point G (2 cm below gluteal fold)

### NB

**c** = circumference, **I** = length  
 All length measurements should be longitudinal and should follow body contours.  
 All measurements to be in centimetres.

Delivery:  to orderer  directly to patient

Name and address of recipient:

Patient Reference \_\_\_\_\_

Name/Last Name \_\_\_\_\_

Clinicians Name \_\_\_\_\_

Street \_\_\_\_\_

Clinicians Work Base \_\_\_\_\_

Postal Code/Place \_\_\_\_\_

Phone \_\_\_\_\_

Phone \_\_\_\_\_

Date of Measurement \_\_\_\_\_

Order Date \_\_\_\_\_

# Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	Extra Large	Extra Extra Large
<b>ReadyWrap® Thigh</b>					
(IG-E1) (2 x lengths)	Groin Circumference (cG) 48-63 cm Above Knee Circumference (cE) 38-48 cm	Groin Circumference (cG) 58-73 cm Above Knee Circumference (cE) 45-58 cm	Groin Circumference (cG) 68-83 cm Above Knee Circumference (cE) 52-68 cm	Groin Circumference (cG) 78-93 cm Above Knee Circumference (cE) 59-78 cm	Groin Circumference (cG) 88-103 cm Above Knee Circumference (cE) 66-88 cm
■ Beige up to 25cm	104478	104474	104470	104482	104790
■ Black up to 25cm	104479	104475	104471	104483	104791
■ Beige up to 30cm	104480	104476	104472	104484	104792
■ Black up to 30cm	104481	104477	104473	104485	104793

<b>ReadyWrap® Knee</b>					
(one length - 20 cm)	Above Knee Circumference (cE) 38-48 cm Below Knee Circumference (cD) 28-38 cm	Above Knee Circumference (cE) 45-58 cm Below Knee Circumference (cD) 33-43 cm	Above Knee Circumference (cE) 52-68 cm Below Knee Circumference (cD) 38-48 cm	Above Knee Circumference (cE) 59-78 cm Below Knee Circumference (cD) 43-58 cm	Above Knee Circumference (cE) 66-88 cm Below Knee Circumference (cD) 53-68 cm
■ Beige	104464	104462	104460	104466	104788
■ Black	104465	104463	104461	104467	104789

<b>ReadyWrap® Calf</b>					
(ID-E) (2 x lengths)	Calf Circumference (cC) 28-38 cm Ankle Circumference (cE) 18-23 cm	Calf Circumference (cC) 33-43 cm Ankle Circumference (cE) 20-28 cm	Calf Circumference (cC) 38-48 cm Ankle Circumference (cE) 25-33 cm	Calf Circumference (cC) 43-58 cm Ankle Circumference (cE) 30-38 cm	Calf Circumference (cC) 53-68 cm Ankle Circumference (cE) 35-43 cm
■ Beige up to 30 cm	104450	104446	104442	104454	104784
■ Black up to 30 cm	104451	104447	104443	104455	104785
■ Beige up to 35 cm	104452	104448	104444	104456	104786
■ Black up to 35 cm	104453	104449	104445	104457	104787

<b>ReadyWrap® Foot</b>					
(Side)	Foot Length (IA-Y) up to 21 cm Ankle Circumference (cE) 18-23 cm Circumference (cY) 28-33 cm	Foot Length (IA-Y) up to 23 cm Ankle Circumference (cE) 20-28 cm Circumference (cY) 30-38 cm	Foot Length (IA-Y) up to 26 cm Ankle Circumference (cE) 25-33 cm Circumference (cY) 35-43 cm	Foot Length (IA-Y) up to 29 cm Ankle Circumference (cE) 30-38 cm Circumference (cY) 40-48 cm	
Right ■ Beige	104436	104434	104432	104438	
Right ■ Black	104437	104435	104433	104439	
Left ■ Beige	104428	104426	104424	104430	
Left ■ Black	104429	104427	104425	104431	

<b>ReadyWrap® Toes</b>					
(cA)	(cA) 18.0-20.4 cm Circumference	(cA) 20.5-23.5 cm Circumference	(cA) 23.6-26.0 cm Circumference	(cA) 26.1-28.5 cm Circumference	(cA) 28.6-31.0 cm Circumference
Right ■ Beige	104416	104414	104412	104418	104420
Right ■ Black	104417	104415	104413	104419	104421
Left ■ Beige	104406	104404	104402	104408	104410
Left ■ Black	104407	104405	104403	104409	104411

<b>ReadyWrap® Extender Strap</b>					
	One size fits all				
■ Beige	104740				
■ Black	104741				

<b>ReadyWrap® Liner extra pairs - thigh length</b>					
(cG) (cC) (cE)		Groin Circumference (cG) 48-83 cm Calf Circumference (cC) 28-48 cm Ankle Circumference (cE) 18-33 cm		Groin Circumference (cG) 78-103 cm Calf Circumference (cC) 43-68 cm Ankle Circumference (cE) 30-43 cm	
■ Black		S - L 149034		XL - XXL 149035	

<b>ReadyWrap® Liner extra pairs - below knee</b>					
(cC) (cE)		Calf Circumference (cC) 28-48 cm Ankle Circumference (cE) 18-33 cm		Calf Circumference (cC) 43-68 cm Ankle Circumference (cE) 30-43 cm	
■ Black		S - L 149032		XL - XXL 149033	

Please mark (tick or cross) where applicable to indicate the garment choice and then enter the sizing, quantity and reference code.

Garment (Left Leg)	Length/Side	Colour	Size	Quantity	Code
Thigh	25 cm 30 cm	Beige Black			
Knee		Beige Black			
Calf	30 cm 35 cm	Beige Black			
Foot	Left	Beige Black			
Toe	Left	Beige Black			
Extender Strap		Beige Black			
Liner		Black			

Garment (Right Leg)	Length/Side	Colour	Size	Quantity	Code
Thigh	25 cm 30 cm	Beige Black			
Knee		Beige Black			
Calf	30 cm 35 cm	Beige Black			
Foot	Right	Beige Black			
Toe	Right	Beige Black			
Extender Strap		Beige Black			
Liner		Black			