


ReadyWrap®

Measuring the lower limb.

How to measure the leg: Ask the patient to stand, if possible. Mark and measure the leg as follows



Write your measurement results in the **box** adjacent to measuring point

Circumference points:

- cG** 2 cm below Gluteal Fold
- cE'** 5 cm above the Knee
- cD** 2-3 cm below Popliteal Fossa
- cC** widest part of Calf
- cB** 2 cm Above Malleolus
- cY** around Instep and Heel with foot dorsi-flexed
- cA** around base of Toes

Lengths:

- IA-Y** heel to base of Toes
- ID-B** the length between point B (2 cm above the malleolus) and point D (2 cm below popliteal fossa)
- IG-E'** the length between point E' (5 cm above the knee) and point G (2 cm below gluteal fold)

NB

c = circumference, **I** = length
 All length measurements should be longitudinal and should follow body contours.
 All measurements to be in centimetres.

Delivery: to orderer directly to patient

Name and address of recipient:

Patient Reference _____

Name/Last Name _____

Clinicians Name _____

Street _____

Clinicians Work Base _____

Postal Code/Place _____

Phone _____

Phone _____

Date of Measurement _____

Order Date _____

Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	Extra Large	Extra Extra Large
ReadyWrap® Thigh					
(G-E1) (2 x lengths)	Groin Circumference (cG) 48-63 cm Above Knee Circumference (cE) 38-48 cm	Groin Circumference (cG) 58-73 cm Above Knee Circumference (cE) 45-58 cm	Groin Circumference (cG) 68-83 cm Above Knee Circumference (cE) 52-68 cm	Groin Circumference (cG) 78-93 cm Above Knee Circumference (cE) 59-78 cm	Groin Circumference (cG) 88-103 cm Above Knee Circumference (cE) 66-88 cm
■ Beige up to 25cm	104478	104474	104470	104482	104790
■ Black up to 25cm	104479	104475	104471	104483	104791
■ Beige up to 30cm	104480	104476	104472	104484	104792
■ Black up to 30cm	104481	104477	104473	104485	104793

ReadyWrap® Knee					
(one length - 20 cm)	Above Knee Circumference (cE) 38-48 cm Below Knee Circumference (cD) 28-38 cm	Above Knee Circumference (cE) 45-58 cm Below Knee Circumference (cD) 33-43 cm	Above Knee Circumference (cE) 52-68 cm Below Knee Circumference (cD) 38-48 cm	Above Knee Circumference (cE) 59-78 cm Below Knee Circumference (cD) 43-58 cm	Above Knee Circumference (cE) 66-88 cm Below Knee Circumference (cD) 53-68 cm
■ Beige	104464	104462	104460	104466	104788
■ Black	104465	104463	104461	104467	104789

ReadyWrap® Calf					
(ID-B) (2 x lengths)	Calf Circumference (cC) 28-38 cm Ankle Circumference (cB) 18-23 cm	Calf Circumference (cC) 33-43 cm Ankle Circumference (cB) 20-28 cm	Calf Circumference (cC) 38-48 cm Ankle Circumference (cB) 25-33 cm	Calf Circumference (cC) 43-58 cm Ankle Circumference (cB) 30-38 cm	Calf Circumference (cC) 53-68 cm Ankle Circumference (cB) 35-43 cm
■ Beige up to 30 cm	104450	104446	104442	104454	104784
■ Black up to 30 cm	104451	104447	104443	104455	104785
■ Beige up to 35 cm	104452	104448	104444	104456	104786
■ Black up to 35 cm	104453	104449	104445	104457	104787

ReadyWrap® Foot					
(Side)	Foot Length (IA-Y) up to 21 cm Ankle Circumference (cB) 18-23 cm Circumference (cY) 28-33 cm	Foot Length (IA-Y) up to 23 cm Ankle Circumference (cB) 20-28 cm Circumference (cY) 30-38 cm	Foot Length (IA-Y) up to 26 cm Ankle Circumference (cB) 25-33 cm Circumference (cY) 35-43 cm	Foot Length (IA-Y) up to 29 cm Ankle Circumference (cB) 30-38 cm Circumference (cY) 40-48 cm	
Right ■ Beige	104436	104434	104432	104438	
Right ■ Black	104437	104435	104433	104439	
Left ■ Beige	104428	104426	104424	104430	
Left ■ Black	104429	104427	104425	104431	

ReadyWrap® Toes					
(cA)	(cA) 18.0-20.4 cm Circumference	(cA) 20.5-23.5 cm Circumference	(cA) 23.6-26.0 cm Circumference	(cA) 26.1-28.5 cm Circumference	(cA) 28.6-31.0 cm Circumference
Right ■ Beige	104416	104414	104412	104418	104420
Right ■ Black	104417	104415	104413	104419	104421
Left ■ Beige	104406	104404	104402	104408	104410
Left ■ Black	104407	104405	104403	104409	104411

ReadyWrap® Extender Strap					
	One size fits all				
■ Beige	104740				
■ Black	104741				

Please mark (tick or cross) where applicable to indicate the garment choice and then enter the sizing, quantity and reference code.

Garment (Left Leg)	Length/ Side	Colour	Size	Quantity	Code
Thigh	25 cm 30 cm	Beige Black			
Knee		Beige Black			
Calf	30 cm 35 cm	Beige Black			
Foot	Left	Beige Black			
Toe	Left	Beige Black			
Extender Strap		Beige Black			

Garment (Right Leg)	Length/ Side	Colour	Size	Quantity	Code
Thigh	25 cm 30 cm	Beige Black			
Knee		Beige Black			
Calf	30 cm 35 cm	Beige Black			
Foot	Right	Beige Black			
Toe	Right	Beige Black			
Extender Strap		Beige Black			