

## measuring & fitting information

**B:** Widest part of the calf

**C:** Just above the malleolus (ankle bone)

**D:** From back of the heel to the longest toe

Activa® British Standard hosiery measuring guide a size indicator can be found by a coloured thread in the stocking					
Size	S	M	L	XL	XXL
<b>B</b>	30.5-37.5cm	33.0-40.0cm	35.5-43.0cm	38.0-46.0cm	40.0-49.0cm
<b>C</b>	19.5-25.5cm	21.5-27.5cm	22.5-29.5cm	23.0-32.0cm	23.5-34.0cm
<b>D</b>	20.5-24.0cm	23.0-26.0cm	25.5-29.5cm	26.5-32.5cm	26.5-32.5cm

